ELx800™ Absorbance Microplate Reader

The ELx800™ is a compact, robust microplate reader ideally suited for applications within the clinical and life science research laboratory. When interfaced with BioTek's Gen5™ software, the ELx800 applications are expanded to include kinetic and well area scanning measurements in microplates from 6-well to 384-well. When space or budget are limited, the ELx800 offers extensive on-board software, complete with multiple curve fit options, data transformations, cutoff and assay validation calculations. The outstanding performance of this hardworking microplate reader and its proven reliability makes the ELx800 an unbeatable value for your laboratory.

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lr@lrdiagnostico.com
www.lrdiagnostico.com

Features:

- Reliable and robust design
- 5 filter capacity
- Multiple plate formats
- Extensive onboard data analysis
- Multiple curve-fitting options
- Assay and control validation
- Data transformation formulas
 - Cutoff and call criteria

Applications:

- ELISA assays
- Enzyme activity

Models:

ELx800: Available for 6- to 384-well microplates and extended dynamic range

See Web site or price list for complete model listings and descriptions.

Optional Accessories:

- Gen5™ Data Analysis Software
- Gen5™CL Custom Clinical Assay Software
- Gen5™ Secure (for 21 CFR Part 11 Compliance)
- Absorbance Test Plate
- Deskjet or Inkjet Printer
- Product Qualification Package

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ELx50™ Microplate Strip Washer

The ELx50™ Microplate Strip Washer's compact footprint conceals a powerhouse of washing capabilities unsurpassed in its class. Exceptional dispense precision and evacuation efficiency can be utilized for both 96- and 384-well strip and plate washing.

Multiple Washers in One

The ELx50 is a 3-in-1 solution for 96-well plate formats automating the wash steps of ELISAs, magnetic bead assays and polystyrene bead assays. The ELx50 can also be equipped with BioTek's Dual-Action™ manifold allowing independent control of dispense and aspiration manifolds

in both 96- and 384-well formats. As a welcome upgrade from manual processing, the ELx50 provides an all-inclusive wash solution offering consistent performance and unattended operation.

Syringe Drive Fluid Delivery

As a self-contained and programmable washer, the ELx50 allows for complete control of precise fluidic delivery from the gentle dripping of a simple squeeze bottle to the full force of pressure delivery systems. Comprehensive onboard software makes creating protocols quick and intuitive.

Automated Liquid Level Sensing

Liquid Level Alert™ allows the convenience of continuous monitoring for both supply and waste bottles. At the beginning and end of a wash protocol, the liquid level is verified to ensure an adequate buffer remains to complete a wash. Sufficient storage capacity in the waste bottle is also verified.

Typical Applications:

- ▶ ELISAs
- ▶ Cell-based assays
- ▶ Magnetic bead assays
- ▶ Polystyrene bead assays
- ▶ ELISPOT assays
- ▶ Multiplex assays

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lr@lrdiagnostico.com

www.lrdiagnostico.com

Specifications

General	
Assays	ELISA Cell-based assays Magnetic bead, polystyrene bead (model dependent) <ul style="list-style-type: none"> • Multiplex assays • Bead-based ELISA Filtration-to-waste processes (model dependent) Protein arrays
Separation	Biomagnetic separation, vacuum filtration (model dependent)
Microplate types	96- and 384-well plates (model dependent) Low profile and standard height Solid and filter bottom (model dependent) <ul style="list-style-type: none"> • Filter pore sizes 0.45 µm to 1.2 µm
Shaking	Programmable in minutes and seconds up to 60 minutes Intensities - 1-5 (15-19 cycles/second)
Soak time	Programmable in minutes and seconds up to 60 minutes
Magnet	High strength 96-well format Flat <ul style="list-style-type: none"> • Flat-bottom well - beads pulled to band across well bottom • Round-bottom well - beads pulled to button at well bottom Ring - beads pulled to 4-zone ring at well bottom
Safety/convenience/maintenance	Waste sensing, fluid detection (optional) Aerosol cover Overflow protection Pre-programmed maintenance routines
Onboard software	Create, edit or run multiple protocols
Washing	
Manifold types	96-well washing: 8-tube (1x8) manifold 12-tube (1x12) manifold 96-/384-well washing: Dual-Action 16-tube (1x16) manifold
Washing speed	Solid bottom plates: 3 asp./disp. cycles: >300 µL/well, 96 wells, 8-tube manifold, final aspirate: ≤130 seconds Filter bottom plates: Variable, based on wash parameters

Fluid delivery	One positive displacement syringe drive
Vacuum filtration	Selectable vacuum levels (final at 30 seconds): 0.45 µm 96-well plates: Low -91 mmHg Medium -150 mmHg High -313 mmHg 1.2 µm 96-well plates: Low -95 mmHg Medium -155 mmHg High -299 mmHg Vacuum filtration time range: 1 - 180 seconds
Volume range	25 - 3,000 µL/well (configuration dependent) Selectable in 1 µL increments
Buffer selection	Automatic switching for up to 3 wash buffers (configuration dependent)
Flow rates	1 - 9
Wash cycles	1 - 10
Dispense precision	≤3% CV
Residual volume	Solid bottom plates: ≤2 µL/well Filter bottom plates: Average increase weight of plate ≤1.2 g after dispensing 300 µL to 0.45 µm plate, filtration 30 seconds, low vacuum, blotted
Sterilization	Chemical: 70% ethyl or isopropyl alcohol and 0.5% sodium hypochlorite solution
Supply bottle	2 L
Physical Characteristics	
Power	100 - 240 Volts AC. 50/60 Hz.
Dimensions	14"W x 16"D x 6.5"H (35.6 x 40.6 x 16.5 cm)
Weight	22 lbs (9.8 kg)
Regulatory	
Regulatory	CE and TUV marked. ROHS compliant. In Vitro Diagnostic use models are available.

Specifications subject to change. Performance values represent the average observed-factory test values.

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BioTek Instruments, Inc.