

New Functions

Display delivers feedback

The eFlow[®]rapid with eBase controller can do much more. Its display delivers feedback during inhalation.




A Few Examples:




Battery charge level 50 %

-> Have new batteries on hand or operate with the power adapter connected




Cable connection to the control unit or on the nebuliser handset interrupted

-> Check cable connections



No medication added

-> Add medication



Good to know:

The eBase controller can be used with the eFlow[®]rapid nebuliser handset for general use as well as with drug-specific nebuliser handsets (e.g. Altera[®] for Cayston[®])



eFlow[®]rapid
nebuliser system

Order No. 178G1005

The eFlow[®]rapid nebuliser handset and the eFlow[®]rapid aerosol head are compatible with both the dark blue Control Unit (eBase Controller) and the turquoise Control Unit (previous version).



eFlow[®]rapid nebuliser handset
(incl. aerosol head)
Order No. 678G8222



eFlow[®]rapid aerosol head
Order No. 678B2620

Local Distributor

The eFlow[®]rapid with eBase Controller



eFlow[®]rapid
nebuliser system

© PARI Pharma GmbH · 178D1051-B · 06/2014



PARI Pharma GmbH
Moosstraße 3 · D-82319 Starnberg

PARI Pharma
Advancing Aerosol Therapies



www.eFlowrapid.info

PARI Pharma
Advancing Aerosol Therapies



**Modern Inhalation Therapy,
it's Twice as Fast**



eFlow[®]rapid enables quick and efficient inhalation significantly reducing the treatment burden^{1,2}.

The Main Advantages at a Glance:

- ➔ Short treatment times mean **better quality of life**
- ➔ Quiet running for **discrete use**
- ➔ Light, compact and portable. Can be run on either batteries or mains power for **greater mobility**
- ➔ Easy to clean, handset can be disinfected and autoclaved for the **highest standards of hygiene**
- ➔ Display **delivers feedback** during inhalation

1) Fischer R et al. 2007, J Cyst Fibr 6 (Suppl. 1): 41.
2) Naehrig S et al. 2011, Eur J Med Res 16: 63-66.

**Hygiene Tips
for home**

Please refer to the Instructions for Use for more details

- ➔ **Clean**
Dismantle the nebuliser into its components. Rinse thoroughly in warm tap water with a little washing-up liquid immediately after every inhalation



- ➔ **Disinfect**
Disinfect the dismantled nebuliser and aerosol head at least once a day by putting them in a commercially-available baby bottle disinfector for 6 minutes or alternatively by boiling them in distilled water for at least 5 minutes



- ➔ **Dry**
Allow to dry completely

- ➔ **Store**
Store in a clean, dry and hygienic place



**Functionality Test
for Aerosol Heads**



Measure the time it takes for 2.5 ml of an isotonic saline solution (NaCl 0,9%) to be nebulised.

- ➔ 1. Should the measured time for nebulisation lie between 4 and 5 minutes, the aerosol head should be **cleaned** using the easycare cleaning aid **4 minutes**

- ➔ 2. Should the measured time for nebulisation be longer than 5 minutes after cleaning with easycare, then the aerosol head should be **changed** **5 minutes**


Components



* Aerosol chamber and medication reservoir are pre-fitted and cannot be taken apart.