


PRODUCT INSTRUCTION SHEET

REFERENCE

Stockerston Blanket Cradle (137)

DESCRIPTION

Blanket support constructed from plastic coated steel tubing. This simple device has closed and open ends, the latter being placed between the mattress and bed base. The weight of the bed clothes is supported, keeping them clean of the patients for ventilation and pressure relief.



TECHNICAL SPECIFICATIONS

| | |
|-------------------------|--------------|
| Overall Width: | 16" (41 cm) |
| Overall Height: | 18" (46 cm) |
| Length of base support: | 30" (76 cm) |
| Product weight: | 1.8 kg (4lb) |

USAGE REQUIREMENTS

N.B No maximum user weight/load is given for this product – this is dependent upon the actual weight of the bed clothes, and we recommend a risk assessment be carried out before use.

1. Ensure the mattress is flat, and the cradle is horizontal or parallel to the top of the mattress.
2. When positioning the unit, make sure it is pushed fully under the mattress and the vertical struts rest against the mattress edge.
3. Check the position of the product regularly to maintain its correct position in line with the original risk assessment
4. Ensure that the coating is not damaged since this could affect hygiene.
5. The product is designed for minimum maintenance, and can be cleaned with a soft cloth using general purpose neutral detergent and warm water. Do not use harsh abrasive materials. To disinfect clean as above and wipe over with a solution of 1000 part per million available chloride disinfectants (haztabs/percept) rinse and dry. If using an autoclave, run at 80C.

Warranty Information

The product warranted for one year from date of purchase due to faulty workmanship and materials. The warranty does not extend to the consequential costs resulting from fault clearance, in particular freight and travel costs, loss of earnings, expenses etc. The manufacturer will not accept responsibility for any damage or injury caused by misuse or non observance of the instructions set out above.