G687004

Gebrauchsanweisung
Operating Instructions
Instrucciones de uso
Manuel d'instructions
Istruzioni d'uso

DE

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EN

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

Medicon implants are manufactured from pure titanium or Titanium alloy. Both materials are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.

- conventional x-ray radiography
- computer tomography
- MRI (Magnetic Resonance Imaging).

The surface is chemically passive and the material is antimagnetic. The implants are supplied in unsterile condition.

Table of contents

1. Indications
 2. Contraindications
 3. Possible side effects and complications
 4. Single Use Product
 5. CT and MR Instructions
 6. Implant selection
 7. Application/handling of plates
 8. Application/handling of screws
 9. Instruments
 10. Decontamination, cleaning and sterilization
 11. Note
- Signs and symbols (See below)

MEDICON eG as the supplier of these products is not liable for damage caused by, or arising from improper use or handling, especially due to non-observance of the operating instructions or improper cleaning or maintenance of the products. The Medicon plate systems may be used only by surgeons with sufficient training and adequate experience in oral, orthodontic and maxillofacial surgery.

1. Indications

The MEDICON Mikro- and Mini implants, consisting of plates and screws, are intended for use in surgical interventions on the osseous structures of the skull. Indications include:

- fractures
- osteotomies, including orthognathic indications
- revision surgery in cases where other methods of treatment or devices were unsuccessful.

The implants have been designed for supporting the normal bone healing process, not for replacing normal body structures or bearing the body weight in cases of insufficient bone healing.

2. CONTRAINDICATIONS

- Patients who are unable to follow the instructions for postoperative care. Such inability may be caused by psychological/mental or neurological conditions of the patient.

- Patients with unstable physical and/or mental health conditions.
- Patients with insufficient or poor-quality bone tissue, circulatory disorders or latent infections.
- Oversensitivity to materials and foreign-body reactions. Preimplantation tests must be carried out even if such oversensitivity or reactions are only suspected.
- Acute infections.

3. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- In many cases, possible complications are attributable not to the implants as such, but to incorrect patient selection, lack of practice or inadequate repositioning of the bone fragments.
- Insufficient adaptation of the implants can lead to bone union delay or failure.
 - Complaints, pains, abnormal sensations or palpability of implants.
 - Material/foreign-body sensitivity of the patient with allergic reactions.
 - Stronger connective-tissue reaction in the fracture region or in the vicinity of the implant.
 - Delayed or insufficient fracture healing leading to possible implant breakage.
 - Inadequate healing.
 - Insufficient bone formation, osteolysis, osteoporosis, reduced revascularization or infection, which can lead to loosening, bending, breakage or rupture of the implant.
 - Breakage, bending, migration or loosening of the implant.
 - Reduced bone density due to "stress shielding".
 - Loosening of the implant as a result of insufficient screw tightening.
 - Osseous necrosis, osteoporosis, inhibited revascularization on, bone-resorption, and poor bone formation, which can cause premature loss of fixation leading to non-union.
 - Extreme or repeated bending of the plates in the same place can result in plate breakage.
 - Unstable comminuted fractures can lead to an increased connective-tissue reaction in the fracture region.
- This surgical procedure may cause not only the above-mentioned side effects and complications but also problems such as injuries to nerves, infections, pain etc., which are not necessarily caused by the implant.

4. SINGLE USE PRODUCT

- Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be reused. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

5. CT & MR INSTRUCTION: MR Conditional

CT- and MR-capable pure titanium plate system are approved for a field strength of 1.5 Tesla when neither warning nor rotation occur. The pa-

tient must be informed by the surgeon of his/her duty to warn the staff about the implant before CT and MRI investigations are performed.

6. IMPLANT SELECTION

The surgeon is responsible for the correct selection of implants.

The following aspects are critical for the correct selection of implants:

- bone defect to be treated / intended bone repositioning,
- body weight of the patient,
- general health condition, degree of activity of the patient.

Incorrect implant selection can lead to premature implant loss and cause loosening, deformations and fractures of the implants. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending, or fracturing of the product and/or bone. The success of an operation also depends on how the implants are handled. The plates must be bent carefully and cautiously, avoiding repeated and excessive bending. Scratches and other damage can reduce the stability of the product and result in early fatigue of the implant. Medicon implants and instruments have been designed for being used together. Using implants or instruments from other manufacturers in conjunction with Medicon products entails unpredictable risks since the products will not be matched to each other. To avoid such risks, only Medicon products that are intended for combination with each other may be used together. Prior to applying the products, the surgeon must thoroughly discuss with the patient the intended outcome of the operation. In such discussions, special attention must be paid to postoperative care measures that may become necessary. The patient must be informed that, due to its limited stability, the fixation must not put under the load of the full body weight. The patient must also be informed that non-observance of this advice can put at risk any safe bone union. The patient must be instructed to inform the surgeon immediately if they notice unusual changes in the situs. The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must consider possible consequences, e.g. implant failure, and discuss with the patient any necessary measures for further healing. In cases of failed, delayed or insufficient fracture union, the incidence of implant bending, implant failure or implant breakage cannot be excluded. Therefore, immobilization of the fracture site must be ensured until there is solid bone fusion. The incessant load changes, to which the implant will be subjected, can lead to fatigue fractures. The implants can be removed once the fracture has healed completely. An implant left in the body after complete healing may act as a load-bearing element and can contribute to an increased risk of re-fracturing in active patients. For this reason, we accept liability for the implants only until healing is complete. We will not accept liability for any damage or injuries caused by implants that were not removed as soon as healing of the treated fracture was complete. The surgeon is responsible for deciding whether the implants shall remain in the body after complete healing of the fracture or whether they shall be removed. For his decision, the surgeon has to weigh risks and advantages of leaving the implants in the body. It cannot be excluded that implants will break, loosen, corrode, migrate through tissue and cause pain. The packaging label carries a "Lot

No.". We recommend copying this "Lot No." to the patient records as this number allows tracing the production history of the implant back to the raw material.

7. APPLICATION AND HANDLING OF PLATES

The plates need to be adapted to the natural or planned bone contours as precisely as possible. Only use the instruments provided by the manufacturer for the purpose of adapting the plates. During the plate bending procedure the plates are subjected to cold stresses, under which the titanium hardens and its ductility is decreased. Therefore it is absolutely necessary to achieve the required plate shape in as few bending steps as possible. Excessive bending can lead to postoperative plate breakage. The plates must not be bent back and forth more than 2 – 3 times. Overaggressive application of the bending instruments can result in visible damage to the implant. When such damage has occurred, a new plate must be selected and adapted more cautiously. Deformed plate holes for the screws entail not only an increased risk of breakage at the holes; they also compromise the precise seating of the screw head in the plate. Therefore the plates must be bent cautiously. Adapted plates must be inspected for notches, deformed screw holes and other mechanical damage before they are applied to the patient. Always select the plates that provide the closest fit to the specific situation and indication. If a plate needs to be shortened, only use the instruments supplied for this purpose. When shortening a plate, care must be taken that the portion of the plate that is cut off is not flung towards the patient or towards third parties during the cutting procedure. The cut edges of a shortened plate must be deburred to prevent tissue injuries or irritations. Bending templates must not be implanted under any circumstances.

8. APPLICATION AND HANDLING OF SCREWS

Self-tapping screws do not require pre-cutting of threads. To apply the drill holes, only use the drills designated by Medicon for the respective screw so that the correct drill hole diameter for the screw is obtained. Only use drill bits with sharp cutting edges! During drilling, there is the risk of heat damage to the bone. Therefore, only work at low drill speed and always ensure sufficient cooling during the drilling procedure. Excessive application of force during the drilling procedure can lead to drill breakage, endangering the patient, the surgeon and third parties. Use the emergency screws intended for this purpose if, due to a failed drilling effort or a worn drill hole, the original screw cannot be anchored firmly in the bone. To ensure that the screwdriver blade is firmly seated in the screw head, axial pressure must be applied when the screwdriver is inserted into the screw head. This ensures the correct axial alignment of screw and screwdriver, and prevents slipping of the screwdriver blade and damage to the screw head. Application of excessive force when tightening the screws can lead to intraoperative screw breakage. Wear of the screwdriver blade will impair the firm connection between blade and screw. Whenever such wear is detected, the screwdriver blade must be replaced by a new blade. Once the implantation has been completed, every screw and plate connection must be checked for stability. If necessary re-tighten the screws. Prior to implant removal, all screw heads must be cleaned carefully, using a scalpel or another suitable instrument, so that the screwdriver can be optimally seated in the screw head.

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