


Each CUI Tissue Expander offers one or more of the following features designed to increase efficiency during use.


Posterior Fill Valve

An exclusive feature on CUI expanders, the Posterior Fill Valve allows the expander to be filled rapidly during surgery. The patented Posterior Fill Valve is available on all CUI Standard, Versafil™, and Specialty expanders with a suggested volume of 50cm³ and greater.



Stainless Steel Connector


Allows adjustments to port tubing length. Supplied with all CUI Standard, Versafil™, and Specialty Expanders.



Remote Injection Ports

Unless otherwise indicated on a specific catalogue item, all CUI Standard, Versafil™ and Specialty Tissue Expanders are supplied with a Low Profile Injection Port. If an alternative design is required, any of the above ports may be exchanged by the surgeon using the Stainless Steel Connector provided. Remote injection ports are designed for use with a 23 gauge or smaller needle. Tubing length is approximately 13cm.


Catalogue Number	Base Diameter (cm)	Projection (cm)
HPP-1	2.9	1.3
LPP-1	2.9	0.8
MDP-1	1.5	0.6



Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

CUI Versafil™ Tissue Expanders

Croissant (8 sizes)



CZV - Croissant Without Stable Base (4 sizes)

Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)	D (cm)
CZV-0907	100	9.0	7.0	4.5	4.5
CZV-1310	300	13.0	10.0	6.5	6.5
CZV-1713	600	17.0	13.0	8.0	8.5
CZV-2116	1100	21.0	16.0	9.5	10.5


CRS - Croissant With Stable Base (4 sizes)

Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)	D (cm)
CRS-0907	100	9.0	7.0	4.5	4.5
CRS-1310	300	13.0	10.0	6.5	6.5
CRS-1713	600	17.0	13.0	8.0	8.5
CRS-2116	1100	21.0	16.0	9.5	10.5

Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

CUI Versafil™ Tissue Expanders

Round (12 sizes)



FZV - Round Without Stable Base (6 sizes)

Catalogue Number	Suggested Volume* (cm³)	A (cm)	B (cm)
FZV-08	200	8.0	6.0
FZV-10	350	10.0	7.5
FZV-12	600	12.0	9.0
FZV-14	1000	14.0	10.5
FZV-16	1500	16.0	11.5
FZV-18	2000	18.0	13.0

SFV - Round With Stable Base (6 sizes)


Catalogue Number	Suggested Volume* (cm³)	A (cm)	B (cm)
SFV-08	200	8.0	6.0
SFV-10	350	10.0	7.5
SFV-12	600	12.0	9.0
SFV-14	1000	14.0	10.5
SFV-16	1500	16.0	11.5
SFV-18	2000	18.0	13.0

* 350cm³ or larger supplied with High Profile Injection port.

Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

CUI Versafil™ Tissue Expanders

Rectangular (14 sizes)



RZV - Rectangular Without Stable Base (7 sizes)

Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)
RZV-0604	60	6.0	4.0	4.0
RZV-0805	125	8.0	5.0	5.0
RZV-1006	225	10.0	6.0	6.0
RZV-1207	350	12.0	7.0	7.0
RZV-1408	500	14.0	8.0	8.0
RZV-1605	350	16.0	5.0	5.5
RZV-2007	650	20.0	7.0	7.7


SRV - Rectangular With Stable Base (7 sizes)

Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)
SRV-0604	60	6.0	4.0	4.0
SRV-0805	125	8.0	5.0	5.0
SRV-1006	225	10.0	6.0	6.0
SRV-1207	350	12.0	7.0	7.0
SRV-1408	500	14.0	8.0	8.0
SRV-1605	350	16.0	5.0	5.5
SRV-2007	650	20.0	7.0	7.7

Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

CUI Speciality Tissue Expanders

SOE - Small Oval (2 sizes)



Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)
SOE-3624	10	3.6	2.4
SOE-4528	20	4.5	2.8


Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

Customer Service

Kronans Droghandel: +46 (0)771 462633


CUI Standard Tissue Expanders

SRS - Rectangular (6 sizes)



Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)
SRS-0604	70	6.0	4.0	3.5
SRS-0805	140	8.0	5.0	4.0
SRS-1006	250	10.0	6.0	5.0
SRS-1206	340	12.0	6.0	5.0
SRS-1207	500	12.0	7.0	6.0
SRS-1508	680	15.0	8.0	6.5

SFS - Round (11 sizes)




Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)
SFS-60	60	6.9	1.4
SFS-80	80	7.6	2.3
SFS-100	100	8.1	3.1
SFS-200	200	10.0	3.9
SFS-300	300	11.4	4.2
SFS-400	400	12.5	5.4
SFS-500	500	13.0	6.3
SFS-600	600	13.5	6.7
SFS-700	700	14.7	6.9
SFS-800	800	15.2	7.1
SFS-1000	1000	16.4	7.6

Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.


CUI Standard Tissue Expanders

SCS - Crescent (3 sizes)



Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)
SCS-0705	80	7.0	4.5	4.0
SCS-1005	250	10.0	8.0	4.0
SCS-1512	500	14.5	10.5	5.0

SLS - Longitudinally Curved (2 sizes)



Catalogue Number	Suggested Volume (cm³)	A (cm)	B (cm)	C (cm)	D (cm)
SLS-1605	300	16.0	9.0	4.0	5.0
SLS-1907	775	19.0	12.0	6.0	7.0

Nominal dimensions: device dimensions are measured with the device placed on a flat surface. Not every device will conform to the dimensions given, slight variations may occur.

CUI

Tissue Expanders

