

## EPILOG PREOP - Supported HD-EEG data formats

EGI/Phillips	*.mff
ANT	*.avr, *.cnt, *.trg
Brain Products	*.eeg, *.edf

## EPILOG PREOP - HD-EEG data requirements

### • Automated Spike Detection

We accept HD-EEG recordings from EGI, Brain Products, ANT or others with a duration > **3 hours** for automated spike detection. If you have multiple files from the same recording, please put the files in a folder first and upload them in a zipped folder.

### • Manual Marking

We accept HD-EEG recordings from EGI, Brain Products, ANT or others with a duration < **3 hours in which the spikes are marked manually.**

The EEG data must be recorded in one **continuous recording session**. We do not except multiple recordings of the same patient in one submission.

If you have multiple files from the same recording, please put the files in a folder first and upload them in a zipped folder.

**Important!** For .mff files of EGI, please zip the files first before uploading.



Questions? contact us: [support\\_preop@epilog.care](mailto:support_preop@epilog.care)