



VHG Labs brand V-Solv™ ICP solvent


Use V-SOLV as a diluent for your oil-based calibration standards and samples for outstanding nebulization characteristics.

V-Solv™ ICP Solvent is our proprietary solvent that is used for diluting oil and other organic liquids for analysis by ICP and ICP-MS. V-Solv™ ICP Solvent offers substantial advantages over conventional low odor/odorless kerosene, xylene and other commercial solvents:

- Very high purity (essentially no trace metals or sulfur) - each bottle is accompanied by a Certificate of Analysis that states the trace impurity levels of 36 metals and sulfur
- Makes very stable dilutions of metallo-organic standards and oil samples
- Extremely low toxicity compared with kerosene or xylene

- Extremely low odour (no kerosene or aromatic smell), which results in a comfortable work environment and is also VOC exempt for consumer products applications (EPA Title 40, Volume 5, Parts 53-59)
- Competitively priced and non-hazardous to ship (flashpoint of 260°F/125°C)
- One gallon bottle size for ease of use and for shipment via common carrier.

Product	Product number
VHG-V-Solv™ ICP Solvent, 1 Gal	V-SOLV-1GAL
VHG-V-Solv™ ICP Solvent, 5 Gallon Cube	V-SOLV-5GAL



For further information, or if you require substances or materials not currently listed please contact one of our local sales offices.

Brazil

Tel: +55 12 3302 5880
Email: bz@lgcstandards.com

Bulgaria

Tel: +359 (0)2 971 4955
Email: bg@lgcstandards.com

China

Tel: 800 810 4630 (toll-free)
Email: infochina@lgcgroup.com

France

Tel: +33 (0)3 88 04 82 82
Email: fr@lgcstandards.com

Germany

Tel: +49 (0)281 9887 0
Email: de@lgcstandards.com

Hungary

Tel: +49 (0)281 9887 0
Email: de@lgcstandards.com

India

Tel: +91 (0)80 6701 2000
Email: in@lgcpromochem.com

Ireland

Tel: +44 7879556983
Email: ire@lgcstandards.com

Italy

Tel: +39 02 2247 6412
Email: it@lgcstandards.com

Middle East

Tel: +49 (0)281 9887 0
Email: de@lgcstandards.com

Netherlands

Tel: +31 (0)643 775 422
Email: nl@lgcstandards.com

Nordic countries

Tel: +49 (0)281 9887 0
Email: de@lgcstandards.com

Poland

Tel: +48 22 751 31 40
Email: pl@lgcstandards.com

Romania

Tel: +40 364 116890
Email: ro@lgcstandards.com

Russia

Tel: +7 812 777 04 88
Email: ru@lgcgroup.com

South Africa

Tel: +27 (0)11 466 4321
Email: info@industrialanalytical.co.za

Spain

Tel: +34 (0)93 308 4181
Email: es@lgcstandards.com

Turkey

Tel: +90 216 360 0870
Email: tur@lgcstandards.com

United Kingdom

Tel: +44 (0)20 8943 8480
Email: uksales@lgcstandards.com

USA + Canada

Tel: 1-855-LGC-USA1 (toll-free)
Email: lgcusa@lgcstandards.com

www.lgcstandards.com • www.vhglabs.com

Science for a safer world

Brazil • Bulgaria • China • France • Germany • Hungary • India • Ireland • Italy • Netherlands
Nordic countries • Poland • Romania • Russia • South Africa • Spain • Turkey • United Kingdom • USA

