



Disc grinder HK200

Mobile but sturdy sample preparation machine for spectrometer samples

- Mobile desk-top model
- For use with self-adhesive grinding paper
- Solid steel cabinet, cast iron disc pan and motor base
- Easy to clean
- Disc pan with exhaust outlet

Technical specifications:

Power supply	230 V / 50 Hz, 370 W
Instrument dimensions (W x D x H)	270 x 290 x 280 mm
Alternative power	230 V / 60 Hz or 115 V / 60 Hz
Mass	22 kg
Revolution speed	2880 rpm
Grinding disc diameter	205 mm

Consumables and accessories

Self adhesive grinding papers

200 mm available as Al_2O_3 , SiC, ZrO_2 in various grades

Sample holders

Magnetic (square)

Magnetic (round)

Mechanical (round)



For more information or to view our reference materials visit www.lgcstandards.com

For further information, or if you require substances or materials not currently listed please contact one of our local sales offices.

Brazil

Tel: +55 12 3302 5880
Email: bz@lgcstandards.com

Bulgaria

Tel: +359 (0)2 971 4955
Email: bg@lgcstandards.com

China

Tel: 800 810 4630
Email: infochina@lgcgroup.com
Territories served:
China, Hong Kong, Macau, Taiwan

Czech Republic

Tel: +420 543 529 205
Email: cz@lgcstandards.com

Finland

Tel: +358 (0)2 233 9355
Email: se@lgcstandards.com

France

Tel: +33 (0)3 88 04 82 82
Email: fr@lgcstandards.com
Territories served:
Algeria, Belgium, France, Luxembourg,
Monaco, Morocco, Tunisia

Germany

Tel: +49 (0)281 9887 0
Email: de@lgcstandards.com

Territories served:

Austria, Cyprus, Germany, Greece, Israel,
Japan, Korea, Singapore, Switzerland,
Thailand, Vietnam

Hungary

Tel: +36 (06) 26 314 891
Email: hu@lgcstandards.com

India

Tel: +91 (0)80 6701 2000
Email: in@lgcpromochem.com

Ireland

Tel: +44 7879556983
Email: ire@lgcstandards.com

Italy

Tel: +39 02 2412 6830
Email: it@lgcstandards.com

Middle East

Tel: +49 (0)281 9887 0
Email: mideast@lgcstandards.com

Netherlands

Tel: +31 (0)643 775 422
Email: nl@lgcstandards.com

Poland

Tel: +48 22 751 31 40
Email: pl@lgcstandards.com

Territories served:

Armenia, Azerbaijan, Belarus, Georgia,
Kazakhstan, Kyrgyzstan, Lithuania, Poland,
Slovakia, Tajikistan,
Turkmenistan, Ukraine, Uzbekistan

Romania

Tel: +40 364 116890
Email: ro@lgcstandards.com

Russia

Tel: +7 (812) 334 48 25
Email: ru@lgcstandards.com

South Africa

Tel: +27 (0)11 466 4321
Email: info@industrialanalytical.co.za

Spain

Tel: +34 (0)93 308 4181
Email: es@lgcstandards.com

Sweden

Tel: +46 (0)33 20 90 60
Email: se@lgcstandards.com

Territories served:

Denmark, Estonia, Iceland, Latvia,
Norway, Sweden

Turkey

Tel: +90 216 360 0870
Email: tur@lgcstandards.com

United Kingdom

Tel: +44 (0)20 8943 8480
Email: uksales@lgcstandards.com
Territories served:
Channel Islands, Ireland, United Kingdom

USA + Canada

Tel: 1-855-LGC-USA1
Email: lgcusa@lgcstandards.com
Territories served:
Canada, Latin America, USA

www.lgcstandards.com