



# Disc grinder HK210

Desktop machine for dry grinding, wet grinding and polishing

- Universal machine for spectrometry and metallography
- Fully adjustable revolution speed with digital display
- Electronically controlled quick-stop brake
- Easy exchange of grinding table with papers of different grades
- Polishing with paper, cloth and abrasive paste, and water as appropriate
- Electromagnetic water inlet valve giving fine regulation of water jet

## Technical specifications:

Power supply	230 V, 50(60) Hz, 400 W
Instrument dimensions (W x D x H)	410 x 280 x 240 mm
Mass	25 kg
Revolution speed	2880 rpm
Grinding disc diameter	200 mm (8 inches)

Revolution speed 280-2000 rpm,  
continuous at constant torque

Water inlet 3/4 inch, max. 6 bar

Water outlet 1 inch

## Consumables and accessories

Spare grinding discs and paper clamping ring

Sample holders

Magnetic and mechanical



For more information or to view our reference materials visit [www.lgcstandards.com](http://www.lgcstandards.com)

**For further information, or if you require substances or materials not currently listed please contact one of our local sales offices.**

### **Brazil**

Tel: +55 12 3302 5880  
Email: [bz@lgcstandards.com](mailto:bz@lgcstandards.com)

### **Bulgaria**

Tel: +359 (0)2 971 4955  
Email: [bg@lgcstandards.com](mailto:bg@lgcstandards.com)

### **China**

Tel: 800 810 4630  
Email: [infochina@lgcgroup.com](mailto:infochina@lgcgroup.com)  
Territories served:  
China, Hong Kong, Macau, Taiwan

### **Czech Republic**

Tel: +420 543 529 205  
Email: [cz@lgcstandards.com](mailto:cz@lgcstandards.com)

### **Finland**

Tel: +358 (0)2 233 9355  
Email: [se@lgcstandards.com](mailto:se@lgcstandards.com)

### **France**

Tel: +33 (0)3 88 04 82 82  
Email: [fr@lgcstandards.com](mailto:fr@lgcstandards.com)  
Territories served:  
Algeria, Belgium, France, Luxembourg,  
Monaco, Morocco, Tunisia

### **Germany**

Tel: +49 (0)281 9887 0  
Email: [de@lgcstandards.com](mailto:de@lgcstandards.com)

#### Territories served:

Austria, Cyprus, Germany, Greece, Israel,  
Japan, Korea, Singapore, Switzerland,  
Thailand, Vietnam

### **Hungary**

Tel: +36 (06) 26 314 891  
Email: [hu@lgcstandards.com](mailto:hu@lgcstandards.com)

### **India**

Tel: +91 (0)80 6701 2000  
Email: [in@lgcpromochem.com](mailto:in@lgcpromochem.com)

### **Ireland**

Tel: +44 7879556983  
Email: [ire@lgcstandards.com](mailto:ire@lgcstandards.com)

### **Italy**

Tel: +39 02 2412 6830  
Email: [it@lgcstandards.com](mailto:it@lgcstandards.com)

### **Middle East**

Tel: +49 (0)281 9887 0  
Email: [mideast@lgcstandards.com](mailto:mideast@lgcstandards.com)

### **Netherlands**

Tel: +31 (0)643 775 422  
Email: [nl@lgcstandards.com](mailto:nl@lgcstandards.com)

### **Poland**

Tel: +48 22 751 31 40  
Email: [pl@lgcstandards.com](mailto:pl@lgcstandards.com)

#### Territories served:

Armenia, Azerbaijan, Belarus, Georgia,  
Kazakhstan, Kyrgyzstan, Lithuania, Poland,  
Slovakia, Tajikistan,  
Turkmenistan, Ukraine, Uzbekistan

### **Romania**

Tel: +40 364 116890  
Email: [ro@lgcstandards.com](mailto:ro@lgcstandards.com)

### **Russia**

Tel: +7 (812) 334 48 25  
Email: [ru@lgcstandards.com](mailto:ru@lgcstandards.com)

### **South Africa**

Tel: +27 (0)11 466 4321  
Email: [info@industrialanalytical.co.za](mailto:info@industrialanalytical.co.za)

### **Spain**

Tel: +34 (0)93 308 4181  
Email: [es@lgcstandards.com](mailto:es@lgcstandards.com)

### **Sweden**

Tel: +46 (0)33 20 90 60  
Email: [se@lgcstandards.com](mailto:se@lgcstandards.com)

#### Territories served:

Denmark, Estonia, Iceland, Latvia,  
Norway, Sweden

### **Turkey**

Tel: +90 216 360 0870  
Email: [tur@lgcstandards.com](mailto:tur@lgcstandards.com)

### **United Kingdom**

Tel: +44 (0)20 8943 8480  
Email: [uksales@lgcstandards.com](mailto:uksales@lgcstandards.com)  
Territories served:  
Channel Islands, Ireland, United Kingdom

### **USA + Canada**

Tel: 1-855-LGC-USA1  
Email: [lgcusa@lgcstandards.com](mailto:lgcusa@lgcstandards.com)  
Territories served:  
Canada, Latin America, USA

[www.lgcstandards.com](http://www.lgcstandards.com)