



PowerStroll Parts List (PWCPP004)

Product Launched: April 2007



The PWCPP004 was launched in April 2007 and has the following characteristics:

- Single drive wheel with black 8" x 1" tyre
- Grey battery bag with 1No. 12V 20ah battery
- Extended plastic push handles (not on all models)
- Right-angled handle
- Lightweight drive unit

Part Code	Description
Z23001-16	Telescopic Tubing Assembly (contains outer tubing, inner tubing and Z23008-16 securement knob)
Z23002-16	Bracketry Assembly (contains pair of brackets for wheelchair)
Z23003-16	Battery Lead for Powerstroll
Z23004-16	Control Box with Lead
Z23006-16	Battery Bag (with Z23003-16 battery lead)
Z23008-16	Securement Knob (for telescopic tubing)
Z23007-16	Hand Control (not shown)
Z23009-16	12V Charger
Z23010-16	PowerStroll Lift Up Strap
Z23011-16	PowerStroll Motor and Leads
Z23012-16	Carry Holdall
Z23013-16	Wheel (with tyre)
Z23014-16	Telescopic Extension Bar (for wheelchairs with 22" seat width)
Z23015-16	Extended Hand Grips (pair)
Z23016-16	Plastic Battery Tray
Z23020-16	Mounting Bracket for 1" (25mm) diameter Push Handle Tubing
PS0014	Shorter Telescopic Tubing Assembly (for wheelchairs with 14" seat width)