


TW003 PARTS MANUAL



ITEM	QTY	NAME	DESCRIPTION
1	1	Z25210-01	TW002/3 RHS frame assembly
2	1	Z25220-01	TW002/3 LHS frame assembly
3	1	Z25230-01	Centre snap link assembly
4	1	Z25240-01	Centre plastic bracket assembly
5	1	Z26310-01	TW002 RHS cable brake assembly
6	1	Z26320-01	TW002 LHS cable brake assembly
7	2	Z25260A01	Rear wheel assembly
8	1	Z25250-01	Front castor assembly
9	1	Z25330-01	TW003 RHS handle assembly
10	1	Z25340-01	TW003 LHS handle assembly
11	1	See Note	Tightening handle assembly

Part code: Z25993
 Applies to 2004 model
 (from batch AE-1817
 [June 2004])

Note: There are two types of tightening handle fitted to the triwalker:
 Z25270-01 Triangular handle with keyed bolt (to clamp tubing)
 Z25170-01 Hexagonal handle with hex head screw (to bolt thru)