

WA002 PARTS MANUAL



ITEM	QTY	NAME	DESCRIPTION
1	1	Z26210-01	WA002 Rear frame assembly
2	1	Z26220-01	WA002 Front frame assembly
3	2	Z26260-01	Rear wheel assembly
4	1	Z26290-01	WA002 LHS handle assembly
5	1	Z26280-01	WA002 RHS handle assembly
6	2	Z26150-01	Front castor assembly
7	2	Z26130-01	Snap link assembly
8	1	Z26140-01	Centre bar assembly
9	1	Z26170-01	Seat assembly
10	1	Z26310-01	RHS cable brake assembly
11	1	Z26320-01	LHS cable brake assembly
12	1	Z26230A01	Foam back rest (new style)
13	1	Z26240-01	Rollator basket
14	1	Z26680-01	Cane holder assembly (not shown)
15	1	Z26182A01	Tightening Handle (new style)

Part code: Z25997
 Applies to 2004 model
 (from batch AE-1817
 [June 2004])