

WA006 PARTS MANUAL



ITEM	QTY	NAME	DESCRIPTION
1	1	Z26610-01	WA006 rear frame assembly
2	1	Z26620-01	WA006 front frame assembly
3	2	Z26640-01	WA006 front castor assembly
4	2	Z26650-01	WA006 rear wheel assembly
5	2	Z26660-01	WA006 cable and brake assembly
6	1	Z26635-01	LHS snap link assembly
7	1	Z26630-01	RHS snap link assembly
8	2	Z26680-01	Tightening knob assembly
9	1	Z25280-01	TW002 RHS handle assembly
10	1	Z25290-01	TW002 LHS handle assembly
11	1	Z26670-01	WA006 PU seat
12	1	Z26675-01	WA006 basket
13	1	Z26680-01	Walking cane holder (not shown)

(Note: Translucent tray only available as an option)

Part code: Z25998
 Applies to 2004 model
 (from batch AE-1817
 [June 2004])