


# WA007 PARTS MANUAL



ITEM	QTY	NAME	DESCRIPTION
1	1	Z26710-01	WA007 Front Frame Assembly
2	1	Z26720-01	WA007 Rear Frame Assembly
3	1	Z26130-01	Snap Link Assembly
4	1	Z26140-01	Centre Bar Assembly
5	2	Z26750-01	WA007 Front Castor Assembly
6	1	Z26760-01	RHS Rear Wheel Ass'y (excludes
7	1	Z26770-01	LHS Rear Wheel Ass'y brake ass'y)
8	1	Z26310-01	RHS Brake Assembly
9	1	Z26320-01	LHS Brake Assembly
10	1	Z26280-01	RHS Handle Assembly
11	1	Z26290-01	LHS Handle Assembly
12	1	Z26170-01	Seat Assembly
13	1	Z26240-01	Basket
14	1	Z26235-01	Curved Backrest
15	2	Z26182A01	Tightening Handle (handles)
16	4	Z26755-01	Tightening Handle (wheels)
17	1	Z26680-01	Cane Holder (two parts)

Part code: Z25995  
Applies to 2005 model