


TITLE:
**WA010 ACTIVE ROLLATOR
 PARTS LIST & DIAGRAM**

DATE: 19/04/05
 DRAWN: PMK
 SCALE: NTS

SHEET No: 1 of 2
 CHECKED:
 REVISION No: —
 DWG No: Z26800

21	Z26821-15	SHOPPING BASKET
20	Z26820-15	BRAKE HANDLE
19	Z26819-15	HAND GRIP
18	Z26818-15	PLASTIC KNOB FOR SPRING, LEFT
17	Z26817-15	PL. KNOB FOR SPR, RIGHT (NOT SHOWN)
16	Z26816-15	SPRING, HANDLE HEIGHT ADJ
15	Z26815-15	FOLDING ROD
14	Z26814-15	HANDLE FOR FOLD MECHANISM
13	Z26813-15	PLASTIC KNOB FOR FOLDING MECH
12	Z26812-15	RED PLASTIC KNOB FOR HANDLE
11	Z26811-15	SPRING, SEAT LOCK
10	Z26810-15	PLASTIC SEAT LOCK, RIGHT
9	Z26809-15	PLASTIC SEAT LOCK, LEFT
8	Z26808-15	HANDLE, SEAT HEIGHT ADJ
7	Z26806-15	SEAT PLATE
6	Z26807-15	HINGE, SEAT PLATE
5	Z26805-15	FORK, FRONT LEFT
4	Z26804-15	FORK, FRONT RIGHT
3	Z26803-15	COVER, BACK WHEEL BOLT
2	Z26802-15	REAR WHEEL
1	Z26801-15	FRONT WHEEL
Item	Name	Description
Parts List – WA010 Active Rollator		

	TITLE:	DATE:	SHEET No:
	WA010 ACTIVE ROLLATOR PARTS LIST & DIAGRAM	19/04/05	2 of 2
		DRAWN: PMK	CHECKED:
		SCALE: NTS	DWG No: Z26800