

Battery Chargers



The Ampage of these chargers ranges from 18A to 15A and can charge batteries from 7Ah to 100Ah.

- » Specifically for deep cycle lead acid batteries
- » 24V output

Compatibility Chart

SCOOTER	BATTERY (SEE NOTE BELOW)		BATTERY CHARGER	
	CODE	SPECIFICATION	CODE	SPECIFICATION
RIO SERIES	12V12AH X2	12V 12Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
KITE	12V12AH X2	12V 12Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
ST1	ELO84BE	12V 12Ah VRLA	51421811800BE	24V D.C. 1.8A UK Plug
SCOUT 12Ah	12V12AH X2	12V 12Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
SCOUT 20Ah	12V17AH X2	12V 17Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
ST1D	ELO84BE	12V 12Ah VRLA	51421811800BE	24V D.C. 1.8A UK Plug
PRISM	12V12AH X2	12V 12Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
PRISM SPORT	12V17AH X2	12V 17Ah VRLA	3D1200300000	24V D.C. 1.8A UK Plug
ST2	ELO84BE	12V12Ah VRLA	51421811800BE	24V D.C. 1.8A UK Plug
ST3	ELO08BE	12V 18Ah VRLA	51422032001BE	24V D.C. 2A UK Plug
NEO	12V34AH X2	12V 34Ah VRLA	C12-001-02100	24V D.C. 5A UK Plug
ENVOY	12V50AH X2	12V 50Ah VRLA	C12-001-02100	24V D.C. 5A UK Plug
ST4E	39618BE	12V 36Ah VRLA	51474058000BE	24V D.C. 5.0A
ST4D	51325073500BE	12V 50Ah VRLA	51474058000BE	24V D.C. 5.0A
REGATTA	12V50AH X2	12V 50Ah VRLA	C12-001-02100	24V D.C. 5A UK Plug
AMBASSADOR	12V50AH X2	12V 50Ah VRLA	C12-001-02100	24V D.C. 5A UK Plug
ST5D	ELO25BE	12V 75Ah VRLA	51428089000BE	24V D.C. 8.0A UK Plug
COBRA	12V70AH X2	12V 79Ah VRLA	C12-005-00800	24V D.C. 8A UK Plug
ST6	ELO25BE	12V 75Ah VRLA	26888BE / 51421592800BE	24V D.C. 15.0A
GT	12V70AH X2	12V 79Ah VRLA	C12-005-00800	24V D.C. 8A UK Plug
KING COBRA	12V100AH X2	12V 96Ah GEL	C12-005-00800	24V D.C. 8A UK Plug
ROYALE	12V70AH X2	12V 79Ah VRLA	C12-005-00800	24V D.C. 8A UK Plug
ROYALE 4 SPORT	12V100AH X2	12V 96Ah GEL	C12-005-00800	24V D.C. 8A UK Plug
SPORT RIDER	12V100AH X2	12V 96Ah GEL	C12-005-00800	24V D.C. 8A UK Plug
EASY RIDER	12V100AH X2	12V 96Ah GEL	C12-005-00800	24V D.C. 8A UK Plug

Note Battery Specification:

Battery specifications are based on the MK battery range from Drive Medical.

They may not be identical to battery already installed on the scooter.

All scooter batteries supplied by Drive Medical comply with IATA Special Provision A67 and are classified as safe for air transportation.