

COBALT

Powerchair

“Highly compact and ultra-manoeuvrable”



CE

The Cobalt Powerchair is a compact, lightweight and highly transportable chair which is simple and easy to use due to its user friendly programmable controller. This powerchair splits into 3 sections for easy handling and storage when not in use.

driveTM
www.drivemedical.co.uk

The Cobalt features a rear wheel drive system which provides the user with optimum manoeuvrability, and along with a tight turning circle, this makes this ultra-responsive powerchair simple and easy to use in any indoor environment.



Splits easily for storage and transportation

Superb Comfort, Style & Performance...

- » Maximum speed of 4mph*.
 - » Maximum range of 8 miles on a full battery charge*.
 - » Available in 46cm (18") seat width.
 - » Designed for indoor and outdoor use.
 - » Compact and manoeuvrable for use in tight spaces and busy environments.
 - » Splits easily for storage and transportation.
 - » Rear wheel drive system for optimum manoeuvrability.
 - » Detachable width adjustable armrests with removable arm pads to suit comfort of user.
 - » Folding backrest for easy storage and transportation.
 - » Height adjustable seat and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences.
 - » Swivel seat facilitates easy transfers on and off the chair.
 - » De-clutch operation to enable you to safely freewheel the chair where and when required.
 - » Padded upholstery for additional comfort and support.
 - » 20cm (8") front castors and 23cm (9") rear wheels with solid puncture proof tyres provide durability and low maintenance.
 - » Height adjustable, flip up footplate to suit users leg length and assist with easy transfers.
 - » Dynamic Shark programmable controller provides an easy and enjoyable driving experience.
 - » Rear anti tip wheels as standard.
 - » Seat belt included for additional safety and security.
 - » Available in two colour options - Metallic Blue or Metallic Red.
 - » Maximum weight capacity 113kg (18st).
- * Subject to conditions of use.



Dynamic Shark programmable controller



Flip up, height adjustable footplate

COBALT

CODE: COBALTX16

Max Weight Capacity:
113kg (18st)



Length:
98cm (38.5")



Height:
33-36cm (13-14")



Width:
61cm (24")



Seat Width:
46cm (18")



Seat Depth:
38cm (15")



Weight:
55kg (121lb)

Colour Options



Metallic Blue

Metallic Red

Leaflet Code: LL211 If you would like to view this product please contact your local stockist, details below.