

SERIAL NUMBER DESCRIPTIONS

COBRA

The serial number sticker for the cobra is located under the seat.

The serial number metal plate is on the left hand side of chassis just in front of the back wheel in line with the batteries.

M48 GT

The serial number sticker is located on the shroud just behind your feet when you sit down.

The serial number metal plate is located just behind the rear left wheel.

Prism Sport (MS019)

The serial number sticker is located on rear right hand side in front of the free wheel lever

The serial number metal plate is located on the back left hand corner of front chassis.

Kite (MS023/24)

The serial number sticker is located underneath the seat

The serial number metal plate is just behind the battery pack.

Envoy (MS050/51/52/58)

The serial number sticker is located underneath the seat

The serial number metal plate is located on the back of the rear chassis.

Scout (MS008)

The serial number sticker is located on the front of the rear chassis

The serial number metal plate is located underneath the scooter on the left hand side (tip)

AVIATOR/ROYALE (HW005)

The serial number sticker is located on the seat post.

There is no serial number metal plate.

Sport Rider

The serial number sticker is located behind your feet when you sit on the scooter

The serial number metal plate is located on the rear chassis (just below the bumper)

Regatta

The serial number sticker is located is on top of the batteries.

The serial number metal plate is located underneath the bumper on the rear chassy

Neo 4/6/8

The serial number sticker is located on the rear shroud underneath the seat

The serial number metal plate is located behind the left hand side rear wheel

Ambassador

The serial number sticker is located just behind the feet on the shroud when you sit down.

The serial number metal plate is located underneath the rear shroud bumper on rear chassy.

Royale

The serial number sticker is locater underneath the ignition.

PCMM08

The serial number sticker is located under the seat on the left hand side behind the mid drive wheel.

The serial number metal plate is located on the front of chassis left hand side behind the footplate.