

CHICAGO

Rise & Recliner /
Standard & Petite




Restwell.

*“Distinct lift,
rest and
recline action.”*

CE

This is without doubt the best value Electric Riser Recliner chair on the market. It displays all the quality associated with the Restwell brand at an excellent price! The Chicago features a distinct lift, rest and recline action including an excellent full snooze position.

www.restwellchairs.co.uk

driveTM
www.drivemedical.co.uk

The Chicago allows you to relax with ease at the simple touch of a button, or alternatively it will gently lift you into a standing position. The Chicago has a stylish design and is upholstered in a Rembrandt fabric or Anti Microbial PVC.

Superb Comfort, Style & Performance...



- » Stylish design to suit any decor.
 - » Ergonomically contoured seat & back and padded upholstery for added comfort and support.
 - » Single motor activates a distinct lift, rest and recline action.
 - » It will rise and tilt users to their feet.
 - » Excellent snooze position with backrest reclining to 64° for supreme comfort and relaxation.
 - » Quiet motors provide a smooth and discreet operation.
 - » Comes with a simple to operate two button handset.
 - » Dual handset model available as a special order (CLR15).
 - » This AM PVC model has been Crib 5 tested, so it is suitable for Nursing, Care and Residential Homes.
 - » Anti-Microbial PVC fabric helps minimise build up of harmful bacteria such as MRSA and is easy to wipe clean for improved infection and hygiene control.
 - » The AM PVC model has a memory foam seat which provides added comfort and helps distribute pressure for long periods in chair.
 - » Quality assured, meeting all relevant fire retardancy standards.
 - » Safety key is built into hand controller to prevent accidental movement and entrapment.
 - » Stitched handset holder to prevent loss or damage to the handset.
 - » Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems.
 - » Available in luxurious Rembrandt fabric with 4 colour options - Terracotta, Gold, Green or Cream.
 - » Available in AM PVC fabric with 2 colour options - Black or Cobblestone.
 - » Maximum weight capacity 150kg (23st).
- * Batteries not included.

CHICAGO

CODE: CLR11 (Standard)

CLR21 (Petite)

CLR15 (AM PVC)

Max Weight Capacity:
150kg (23st)




 **Height:**
108cm (42½") - ST
100cm (39½") - PE

 **Width:**
83cm (32¾") - ST
77cm (30¼") - PE

 **Depth:**
93cm (36½") - ST
92cm (36") - PE

 **Seat Depth:**
52cm (20½") - ST
48cm (19") - PE

 **Seat Height:**
48cm (19") - ST
44cm (17½") - PE

 **Seat Width:**
50cm (20") - ST
47cm (18½") - PE

Rembrandt Fabric (Standard & Petite)



Green



Gold



Terracotta



Cream

Anti Microbial PVC (Standard Only)



Cobblestone



Black

Leaflet Code: LL126

If you would like to view this product please contact your local stockist, details below.

