

# INTEGRA

Shell Seat For Nursing Homes

*“Designed to  
provide exceptional  
manoeuvrability”*



CE

Developed as an aid for moving and handling residents around in nursing home environments. It is an ideal transfer aid and provides the user with great comfort and support.

**drive**<sup>TM</sup>  
[www.drivemedical.co.uk](http://www.drivemedical.co.uk)

The Integra Shell Seat tilts forward easily to help the user stand up from the chair. The chair is useful for users who are seated for long periods at a time and at risk of developing pressure ulcers, therefore requiring additional support, comfort and positioning assistance.



## Superb Comfort, Style & Performance...

- » Tilt in space chair - gas strut operated.
- » The chair will rise and tilt users to their feet.
- » Easy inclination due to the ergonomically designed push handle with hand lever.
- » Angle adjustable and removable legrest which automatically moves under the seat when the chair is tilted forwards, enabling the user to get out of the chair with ease.
- » Velcro attached ergonomically designed headrest, which can be easily adjusted.
- » Lap belt is fitted as standard for additional safety.
- » Handy storage pocket on the backrest.
- » Easy to manoeuvre high quality brake wheels (12.5cm) which integrally rotate 360 degrees.
- » Directional wheel enables the chair to be transferred in a straight line.
- » Good hoist access for easy moving and handling.
- » Brown leather style upholstery on the outside, together with 2 way stretch, breathable Dartex material on the inside.
- » High risk cushion supplied as standard.

CODE: SHELLINT16  
(Integra Shell Chair, 16" / 39cm)

CODE: SHELLINT18  
(Integra Shell Chair, 18" / 45cm)

CODE: SHELLINT20  
(Integra Shell Chair, 20" / 51cm)

Max Weight Capacity:  
120kg (19st)



Seat Width:  
39, 45, 51cm  
(15, 18, 20")



Width:  
67 - 72cm  
(26 - 28")



Seat Depth:  
41, 44, 47cm  
(16, 17, 18.5")



Armrest Height:  
23cm (9")



Seat Height:  
45, 49, 53cm  
(18, 19, 21")



Shell Angle:  
-3° Transfer Forward  
+27° Backwards



Depth:  
89cm (35")



Seat Weight:  
from 29 kg (4.5st)



Ergonomic wing headrest



Angle adjustable legrest



Easy inclination thanks to the push-handle shape and the adjustable trigger position.



Sewn storage bag at the back

Leaflet Code: LL215 If you would like to view this product please contact your local stockist, details below.

