



Silver Sport Wheelchair - Spare Parts Catalogue

Released 2009: Revised 2010

2010 Version listed, contact Drive Medical for earlier models.

A simple and reliable wheelchair with attractive silver vein finish and plastic mag wheels with solid tyres.

Comes with removable footplates and padded armrests

CONSUMABLES AVAILABLE

Z11000-24	FOOTREST END CAP, SSP118FA-SF
Z11001-24	CASTOR HOUSING CAP, SSP118FA-SF
Z11002-24	RUBBER BRAKE HANDLE, SSP118FA-SF

Silver Sport Wheelchair (SSP118FA-SF) Spares List

<i>Qty</i>	<i>Code</i>	<i>Description</i>
1	Z11010-24	Centre Pin Piece
1	Z11020-24SV	SSP118FA-SF RHS Frame Assembly
1	Z11030-24 SV	SSP118FA-SF LHS Frame Assembly
1	Z1113 0-24SV	SSP118FA-SF RHS Cross Brace Assembly
1	Z1114 0-24SV	SSP118FA-SF LHS Cross Brace Assembly
1	Z11210-24	18" Seat Canvas Assembly
1	Z112 60-24	18" Back Canvas Assembly
2	Z11400-24	Castor Wheel and Fork Assembly
2	Z11450-24	Nut and Bolt for Rear Axle
2	Z11460-24MAG	24" Rear Mag Wheel Assembly
2	Z10592-24	Economy arm pad
1	Z11600-24BLK	SSP118FA-SF RHS Leg Rest Assembly
1	Z11650-24 BLK	SSP118FA-SF LHS Leg Rest Assembly
1	Z11700-24	SSP118FA-SF RHS Brake Assembly
1	Z11750-24	SSP118FA-SF LHS Brake Assembly
2	Z11024-24	Plastic Seat Guide - Front
2	Z1102 5-24	Plastic Seat Guide - Rear

Note: For handed parts, RHS is the right hand side of the wheelchair when sat in the chair and LHS is the left hand side when sat in the chair.

drive
MEDICAL DESIGN & MANUFACTURING

Ainsley's Industrial Estate, Eiland HX5 9JP
T: +44 (0) 1422 31 4488
E: enquiries@drivemedical.co.uk
W: www.drivemedical.co.uk/technical

TITLE:

SSP118FA-SF
SILVER SPORT WHEELCHAIR
PARTS CATALOGUE

DATE:

29-06-09

DRAWN:

PMK

CHECKED:

SCALE:

NTS

SHEET No:

REVISION No:

DWG No:

SSP118FA