

*This parts catalogue shows the electrical parts available for Restwell Rise Recline and Recline only armchairs. It shows the parts available for each model, with a quick reference list and part photos for illustration:*



**For Boston Rise Recline Armchairs – Models with Jacquard Weave Fabric**

**(Model code: CLR3Jx, CLR1Jx (Petite) )**

Applies to models made with jacquard weave fabric (with diamond pattern) which were made up to May 2013.

- **PACKA** (includes handset, actuator, extension cable and control box)
- **Z50829-13**: Standard **TA1** Actuator 210mm
- **Z50885-13**: Standard **TH1** Two Button Handset
- **Z50850-13**: Standard **TC2** Control Box
- **Z50881-13**: Key Jack and Chain
- **Z50810-13**: Extension Cable for Handset
- **FRAMEPACK3** (Contains LH and RH chair scissor mechanisms)



**For Boston Rise Recline Armchairs – Models with Sophia Fabric**  
**(Model code: CLR3Yx, CLR1Yx – Petite) – Direct out motor set**

Applies to models made with dobby Sophia fabric (with dotted pattern) which were made after June 2013.

- **PACKS** (includes **Economy 2-button** handset, **Contain** actuator and **power supply with battery backup** )
- **Z50838-13**: Economy **TA6** Actuator 210mm
- **Z50897-13**: Economy **TH6** 2-button Handset
- **Z50858-13**: Power Supply **TP2 with Battery Backup** and Mains Cable
- **Z50881-13**: Key Jack and Chain
- **FRAMEPACK3** (Contains LH and RH chair scissor mechanisms)



For New York Snooze Armchairs (Model code: CLR5x)

- **PACKA** (includes handset, actuator, extension cable and control box)
- **Z50829-13: Standard TA1 Actuator 210mm**
- **Z50885-13: Standard TH1 Two Button Handset**
- **Z50850-13: Standard TC2 Control Box**
- **Z50881-13: Key Jack and Chain**
- **Z50810-13: Extension Cable for Handset**



For Denver Dual Motor Armchair (Model code: CDL9x)

- **PACKB** (Includes 2x actuators, 5-button handset and control box)
- **Z50825-13: Standard TA1 Backrest Actuator 172mm**
- **Z50820-13: Standard TA1 Legrest Actuator 185mm**
- **Z50880-13: Standard TH1 5-button Handset**
- **Z50850-13: Standard TC2 Control Box**
- **Z50881-13: Key Jack and Chain**
- **Z50810-13: Extension Cable for Handset**

- **FRAMEPACK4** (Contains LH and RH Scissor Mechanism)

For Frisco Spacesaver Chair (Model Code: CLR6x)



- **PACKC** (includes handset, actuator, extension cable and control box)
- **Z50828-13: Standard TA1 Frisco Actuator 203mm**
- **Z50885-13: Standard TH1 Two Button Handset**
- **Z50850-13: Standard TC2 Control Box**
- **Z50881-13: Key Jack and Chain**
- **Z50810-13: Extension Cable for Handset**



For Chicago Riser Chair In Fabric and PVC (Model Code: CLR11x)  
For Chicago Petite Riser Chair In Fabric and PVC (Codes: CL21x)

- PACKQ (Contains actuator, Power supply, economy 2-button handset)
- Z50838-13: Economy TA6 Actuator 210mm
- Z50897-13: Economy TH6 2-button Handset
- Z50856-13: Power Supply TP5 for Single Motor and Mains Cable
- Z50881-13: Key Jack and Chain
  
- FRAMEPACK1 (Contains LH and RH chair scissor mechanisms)

*Note: older chairs may be fitted with PACKD or PACKO electrics.*



For Chicago Riser Chair In PVC with two handsets (Model Code: CLR15x)

- PACKJ (Contains actuator, control box, 2-button handset and 3-button handset)
- Z50834-13: Economy TA6 Actuator 210mm
- Z50890-13: Economy TH6 2-button Handset
- Z50887-13: Standard TH13-button Handset
- Z50853-13: Control Box TC2B for Single Motor
- Z50881-13: Key Jack and Chain
  
- FRAMEPACK1 (Contains LH and RH chair scissor mechanisms)



Nevada and Portland Dual Motor Armchairs (Model Code: CLR12x)  
Nevada and Portland Petite Dual Motor Armchairs (Code: CLR22x)

- PACKF (Contains 2x actuators, control box, and economy four button handset)
- Z51270-13: Economy TA6 Leg Actuator 210mm
- Z51275-13: Economy TA6 Backrest Actuator 132mm
- Z51280-13: Economy TH6 4-button Handset
- Z51250-13: Control TC2B Box for Dual Motor
- Z50881-13: Key Jack and Chain
  
- FRAMEPACK2A (Contains LH and RH footrest scissor)

mechanisms – CLR12 Standard Size)

- FRAMEPACK22A (Contains LH and RH footrest scissor mechanisms – CLR22 Petite Size)



**For Houston Button Back Chair (Model Code: CLR14x)**

- PACKQ (Contains actuator, Power supply, economy 2-button handset)
- Z50838-13: Economy TA6 Actuator 210mm
- Z50897-13: Economy TH6 2-button Handset
- Z50856-13: Power Supply TP5 for Single Motor and Mains Cable
- Z50881-13: Key Jack and Chain
  
- FRAMEPACK1 (Contains LH and RH chair scissor mechanisms)

*Note: older chairs may be fitted with PACKD or PACKO electrics.*



**For Dallas Button Back Chair (Model Code: CLR18x)**

- PACKQ (Contains actuator, Power supply, economy 2-button handset)
- Z50838-13: Economy TA6 Actuator 210mm
- Z50897-13: Economy TH6 2-button Handset
- Z50856-13: Power Supply TP5 for Single Motor and Mains Cable
- Z50881-13: Key Jack and Chain
  
- FRAMEPACK1 (Contains LH and RH chair scissor mechanisms)



**For Indiana Standard and Petite Chairs (Model Code: CLR19 and CLR29x)**

- PACKQ (Contains actuator, Power supply, economy 2-button handset)
- Z50838-13: Economy TA6 Actuator 210mm
- Z50897-13: Economy TH6 2-button Handset
- Z50856-13: Power Supply TP5 for Single Motor and Mains Cable



- Z50881-13: Key Jack and Chain
- FRAMEPACK1 (Contains LH and RH chair scissor mechanisms)



For Oakland Chair with Knuckles (Model Code: CLR23x)

- PACKS (Contains actuator, power supply with battery backup, economy 2-button handset)
- Z50838-13: Economy TA6 Actuator 210mm
- Z50897-13: Economy TH6 2-button Handset
- Z50858-13: Power Supply TP2 with Battery Backup and Mains Cable
- Z50881-13: Key Jack and Chain



Memphis Dual Motor Armchair (Model Code: CLR24x / CLR24x)

- PACKF (Contains 2x actuators, control box, and economy four button handset)
- Z51270-13: Economy TA6 Leg Actuator 210mm
- Z51275-13: Economy TA6 Backrest Actuator 132mm
- Z51280-13: Economy TH6 4-button Handset
- Z51250-13: Control TC2B Box for Dual Motor
- Z50881-13: Key Jack and Chain
- FRAMEPACK2A (Contains LH and RH footrest scissor mechanisms - CLR24 Standard Size)
- FRAMEPACK22A (Contains LH and RH footrest scissor mechanisms - CLR34 Petite Size)



For Intalift Armchairs (Model Code: ILR4x, ILR5x, ILR6x)

*Aug 2008 - Dec2014 (for models manufactured before Aug 2008 contact Drive Medical)*

- PACKG (Contains 3x actuators, 5-button handset, control box and power supply)
- Z50894-20: TA6 Backrest Actuator 89mm
- Z50893-20: TA6 Legrest Actuator 240mm
- Z50942-20: TA6 Lift Actuator 225mm
- Z50896-20: TH1 5-button Handset without Key Jack
- Z50891-20: TC3 Control Box without Power Supply



- **Z50895-20: TP3 Power Supply with Battery Back Up For Seattle (ILR4) and Dakota (ILR12) Intalift Chairs 2014 onwards**

**PACKZ** (Contains 3x actuators, 5-button handset, control box and power supply) :

- **Z50894-13: TA6 Backrest Actuator 89mm (ILR4 and ILR12)**
- **Z50893-13: TA6 Legrest Actuator 240mm (ILR4 and ILR12)**
- **Z50942-13: TA6 Lift Actuator 225mm (ILR4 and ILR12)**
- **Z50896-13: TH1 5-button Handset without Key Jack (ILR4 and ILR12)**
- **Z50891-13: TC3 Control Box without Power Supply (ILR4 and ILR12)**
- **Z50895-13: TP3 Power Supply with Battery Back Up (ILR4 and ILR12)**



#### Washington Chair

##### Single Motor Riser Chair (CLC001)

**PACKAA** (Containing Actuator, 2-button Handset and Power Supply):

**Z50823-51: Actuator**

**Z50150-51: Power Supply (with mains lead and battery back up)**

**Z50162-51: 2-button Handset with Key Jack**

##### Dual Motor Riser Chair (CLC002)

**PACKAB** (Containing 2x Actuators, 6-button Handset and Power Supply)

**Z50821-51: Legrest Actuator**


**Z50822-51: Backrest Actuator**

**Z50150-51: Power Supply (with mains lead and battery back up)**

**Z50161-51: 6-button Handset with Key Jack**



*Lars pictured above, Sven below*



### Lars and Sven Leather Armchairs

#### Recline Only Chairs (ALC2105x, ALC2117x)

- PACKH (includes handset, actuator, extension cable and control box)
- Z50831-21: TA5 Actuator
- Z50886-21: TH6 2-button Handset (with key jack)
- Z50853-21: TC2B Control box

#### Rise Recline Chairs – Single Motor (ALC2205x, ALC2217x)

- PACKI (includes handset, actuator, extension cable and control box)
- Z50830-21: TA6 Actuator
- Z50885-21: TH1 2-button Handset (with key jack)
- Z50851-21: TC2B Control box
- Z50805-01\* Lars backrest (please state colour)

#### Rise Recline Chairs – Lars Dual Motor (ALCLD2205x)

- PACKM (includes handset, actuators, extension cable and control box)
- Z50834-21: TA6 Backrest Actuator
- Z50835-21: TA6 Legrest Actuator
- Z50890-21: TH1 4-button Handset (with key jack)
- Z50855-21: TC2B Control box

#### Rise Recline Chairs – Borg Single Motor (BRC001)

- PACKX (includes handset, actuator, and control box)
- Z50842-49: TA6 Backrest Actuator
- Z50897-49: 2-button Handset (with key jack)
- Z50865-49: Power Supply with Battery Backup

#### Rise Recline Chairs – Borg Dual Motor (BRC002)

- PACKY (includes handset, actuators, and control box)
- Z50871-49: Four button handset
- Z50843-49: TA6 actuator – leg rest
- Z50844-49: TA6 actuator – backrest
- Z50865-49: Power Supply with Battery Backup
- Z50859-49: Control Box (TJ2-1151-001)



Toronto PVC Chairs

Recline Only Armchairs (Model code: GCERx)

- PACKK (Contains actuator, 2-button handset and control box)
- Z50833-21: Actuator
- Z50889-21: 2-button Handset
- Z50854-21: Control Box

Rise Recline Armchairs (Model code: GCRRx)

- PACKL (Contains actuator, 2-button handset and control box)
- Z50832-21: Actuator
- Z50888-21: 2-button Handset
- Z50852-21: Control Box