


For chairs
manufactured from
01/09/2009 c/w T.
Motion Actuators.

Anti A Entrapment
supplied on chairs
from October 2010



Item desc	Code	Spec desc	Image	Applicable
actuator	ZH920101BE	for OEC5 , stroke 210mm , assy size 388mm , push 6000N , Pitch 3.5 , support AED , ext lead 200mm		Serena
actuator	ZH922101BE	for WM frame , stroke 192mm , assy size 373mm , push 6000N , Pitch 3.5 , support AED , ext elad 200mm		Sofia/Sienna
actuator	ZH923101BE	for OEC2 backrest , stroke 183mm , assy size 371mm , push 6000N , support AED , ext lead 500mm		Serena dual motor
actuator	ZH923102BE	for OEC2 seat , stroke 184mm , assy size 410mm , push 6000N , support AED , ext lead 500mm		Serena dual motor
Cable hub (ECU - Electronic Control Unit)	ZH920201BE	Common cable hub for single/dual motor, support AED, DC cable 2000mm, not include mounting kits		Serena/Sofia/Serena dual motor/Sienna
Mounting kit	ZH923201BE	mounting cable hub to the frame for dual motor chairs		Serena dual motor

Item desc	Code	Spec desc	Image	Applicable
Inline power box	ZH920301BE ZH920302BE ZH920303BE	AC230V , DC cable 100mm , AC cable 1500 detachable		Serena/Sofia/Serena dual motor and Sienna
SMPS power box	ZH920311BE ZH920312BE ZH920313BE	AC100V~240V , DC cable 100mm , AC cable 1500 detachable		Serena/Sofia/Serena dual motor/Sienna
handset dual motor	ZH923401BE	Dual motor handset 5 buttons , support AED , cable 900mm(flexi 300mm) , QR 180° connector, print chair image on handset		Serena dual motor
handset single motor	ZH920401BE	Single motor handset 2 buttons , support AED , cable 900mm(flexi 300mm) , QR180°connector, print chair image on handset		Serena/Sofia/ Sienna
handset ext lead	ZH920402BE	handset extension lead , QR , cable 1400mm		Serena/Sofia/Serena dual motor/ Sienna
Anti entrapment device 920, 922, 924	ZH920501BE	3 AL sections 58x50x58cm. Total 2000mm (without leads) cable 1000mm, connecting leads 300mm, with assembling brackets from T motion, screws, cable fixing blocks and cable ties from ZH		Serena/Sofia/Sienna
Anti entrapment device 923	ZH920502BE	3 AL sections 72x54x72cm. Total 2000mm (without leads) cable 2000mm, connecting leads 300mm, with assembling brackets from T motion, screws, cable fixing		Serena Dual Motor
Anti Entrapment c/w fixings for Sofia	AE100	Anti Entrapment c/w fixings for Sofia		Sofia

Item desc	Code	Spec desc	Image	Applicable
Anti Entrapment c/w fixings for Serena	AE200	Anti Entrapment c/w fixings for Serena		Serena
Anti Entrapment c/w fixings for Serena Dual	AE300	Anti Entrapment c/w fixings for Serena Dual		Serena dual motor
Anti Entrapment c/w fixings for Sienna	AE400	Anti Entrapment c/w fixings for Sienna		Sienna
Anti Entrap. Curtain	ZH920503BE	With velcro and straps	n/a	Serena
Anti Entrap. Curtain	ZH920504BE	With velcro and straps	n/a	Sofia, Sienna
Anti Entrap. Curtain	ZH920505BE	With velcro and straps	n/a	Serena Deluxe
Anti Entrap. Curtain	ZH925505BE	With velcro and straps	n/a	Sasha
923 Castor set	EL0112BE	Diameter 40 x 40		923 Dual Motor - common
924 back castor set	EL0114BE	M8 x 15, 45mm		Sienna - rear
924 front castor set	EL0115BE	58mm		Sienna - front