


# Transfer Turntable

- Easy turning movement.
- Helps nursing person work easier.
- Prevent back strains.



This strong and lightweight turntable provides easy rotation round 360°. Use with a person who can stand but has difficulty turning when making a transfer. It helps nursing person the safe transport of person from chair to bed or to car etc. and helps to avoid back strains.

Allowable carrying load: 115kgs, weight: 1kg, 405mm in diameter.

**DO NOT SIT ON THIS TURNTABLE, IT IS ONLY DESIGNED TO STAND ON.**



CE  
Made in Taiwan

