



No.	Code	Description
1	C10-046-00200	VR2 Power Module
2	C10-028-00100	VR2 Joystick
3	S03-028-00101	Front Fork Assy, Geo
4	S45-028-00101	Footrest Assy, Geo
5	C07-028-00300	Seat Post, Geo
6	C06-028-00100	Rear Wheel with Grey Tyre, Geo
7	S08-028-00500	Front Shroud Assy, Geo
8	S08-028-00618	Rear Shroud Assy, Red, Geo
8	S08-028-00654	Rear Shroud Assy, Blue, Geo

No.	Code	Description
9	C05-028-00100	Front Wheel w. Grey Tyre, Geo
10	C07-028-00404	Armrest, Left, Geo
11	S07-028-00500	Armrest, Right, VR2, Geo
12	RIOLBP	12ah Pack W/O Batteries
13	S09-028-00201	Motor Assy, Right, Geo
14	S09-028-00101	Motor Assy, Left, Geo
15	S07-028-00100	Seat, Beige, Geo
	C07-028-01000	Lap Strap

\*<sub>1</sub> Note: Battery Pack includes batteries

\*<sub>2</sub> The Geo (S3500) powerchair was previously called the Micro Mambo 211 Powerchair (PCMM04).