

MOUNTWAY SOLO®

Toilet Lift

“Powered lifting
that replicates the
body’s natural
movement”



CE

The Mountway Solo® Toilet Lift replicates natural standing movement at the touch of a button, assisting users to help themselves get up into a standing position.

www.drivedevilbiss.co.uk

drive
DeVilbiss®
HEALTHCARE

The Mountway Solo® Toilet Lift is ideal if you struggle to get on or off the toilet, perhaps due to muscle weakness or stiff joints. The lift will lower you gently onto the toilet, then bring you effortlessly back up into a natural standing position.



Powered lifting action replicates body's natural standing movement



Fits easily and compactly over existing toilet, blending discreetly into any home



Arms with or without an integral operating switch offer additional support

Features and Benefits

- › Battery powered device that enables users with restricted mobility to get on and off toilet at the touch of a button
- › Ideal for individuals who have difficulty performing sit to stand activities due to muscle or joint weakness in their legs and/or upper body e.g. due to multiple sclerosis, osteo or rheumatoid arthritis, hip replacement, heart problems and obesity
- › Helps reduce manual handling risks for carers who provide toileting assistance
- › Powered lifting action replicates body's natural standing movement
- › Seat can be stopped at any position to suit user comfort and support requirements
- › Fits easily and compactly over existing toilet, blending discreetly into modern home
- › Arms with or without an integral operating switch offer additional support. They tilt up for easy side transfers or for others to us the toilet
- › Smooth, wipe clean surfaces for ease of cleaning
- › Acts as a normal toilet seat in lowered position for use by other household members
- › No specialist installation - simply remove existing toilet seat and place over toilet bowl
- › Feet rotate clockwise and anti clockwise to adjust height for different users/toilets
- › Models also available with tilt up armrests (with/without integral operating controls), which give additional user support
- › Supplied with two rechargeable batteries to ensure continuous operation
- › Battery level indicator monitors when recharging is required
- › Set of extended feet is available as an accessory to raise seat a further 6cm (2.5") for accommodating tall users or high toilets (Product code: TL634)

Specifications:

CODES: MS0HSUK (Solo Toilet Lift without Arms)
 MS2LSUK (Solo Toilet Lift with 2 Arms (Left Switch))
 MS2RSUK (Solo Toilet Lift with 2 Arms (Right Switch))

Width

61cm (24") (w/o Arms)
 71cm (28") (with Arms)

Depth:

54cm (21")

Height:

43cm – 49cm (17" – 19") (w/o Arms)
 66cm – 72cm (26" – 28") (with Arms)

Seat Width:

55cm (22")

Seat Down Height:

43cm – 49cm (17" – 19")

Seat Up Height:

56cm – 62cm (22" – 24")

Pedestal Clearance Height:

36cm – 42cm (14" – 17")

Toilet Clearance Width:

39cm (15")

Weight (Excl. Battery):

18kg (40lb) (w/o Arms)
 20kg (45lb) (with Arms)

Weight Capacity:

160kg (25st)

Leaflet Code: LL238

If you would like to view this product please contact your local stockist, details across.