



1. FRAME, TITAN18LC  
CODE: C01-072-00100
2. SEAT POST, TITAN18LC  
CODE: C07-072-00100
3. FOOTREST MOUNTING BRACKET, TITAN18LC  
CODE: C45-072-00201
4. FOOTREST, TITAN18LC  
CODE: C45-072-00300
5. PU FOAMED TIRE, 260 X 85, BLACK  
CODE: C06-053-00301
6. CASTOR WHEEL, TITAN  
CODE: C05-072-00100V
7. BATTERY CASE, TITAN18LC  
CODE: C08-072-00500
8. SHROUD, TITAN18LC  
CODE: C08-072-00100
9. 12VOLT 34AMP AGM BATTERY 2 PER OUTER CARTON MU-1 SLD M-2  
CODE: 12V34AH
10. ANTI-TIP WHEELS  
CODE: C28-032-00100
11. MOTOR ASSY (CW MOTOR, EM BRAKE, GEARBOX), RIGHT, TITAN18LC  
CODE: C09-072-00100
12. MOTOR ASSY (CW MOTOR, EM BRAKE, GEARBOX), LEFT, TITAN18LC  
CODE: C09-072-00200
13. CONTROLLER SHARK 40A, UPPER (DK-REMD01)  
CODE: C10-072-00100
14. CONTROLLER SHARK 40A, LOWER (DK-PMB40)  
CODE: C10-072-00200
15. TITAN ARMREST, ADJUSTABLE, L  
CODE: C07-072-00200
16. TITAN ARMREST, ADJUSTABLE, R  
CODE: C07-072-00300
17. ARM PAD, TITAN  
CODE: C07-072-00400
18. SEAT, LOW BACK FOR TITAN18LC  
CODE: C07-072-00700
19. MAIN WIRE HARNESS, TITAN18LC  
CODE: C18-072-00100
20. BATTERY CABLE (+), TITAN18LC  
CODE: C18-072-00300
21. BATTERY CABLE (-), TITAN18LC  
CODE: C18-072-00200
22. DECORATION SET, SPIRIT RED 18, TITAN18LC  
CODE: S08-072-00118