





2016

Eclipse II & Elite Parts List


Eclipse II (PHP398)	1006060001
Elite (PHP399) ▲	1006040001




- | | | |
|----|--|------------|
| 1 | KNOB Comfort Control Eclipse II | 1009060006 |
| 2 | PANEL Facia Complete for Eclipse II Pump | 1009060013 |
| 2a | PANEL Control Front only Eclipse II | 1009060011 |
| 3▲ | LABEL Weight settings Comfort dial (Green) Elite | 6004990040 |
| 4 | CASE Front Blank for Eclipse II Pump | 1009060004 |




- | | | |
|---|--|------------|
| 5 | CASE Back for Eclipse II Pump | 1009060003 |
| 6 | BEDHOOK inc Cover PVC Black Prelude / Eclipse I/II / Elite | 1009980035 |
| 7 | BRACKET Bedhook Prelude / Eclipse I/II / Elite | 1009980036 |
| 8 | SCREW Bedhook M3x8mm Prelude / Eclipse I/II / Elite | 1009980034 |
| 9 | FOAM Black Strip Anti Rattle Pump 1000 L18cmxW1.5cmxH0.8cm | 6003990084 |




- | | | |
|----|--|------------|
| 10 | SCREW Back Case for Phase3/Pump1000/Eclipse | 1009980094 |
| 11 | BEDHOOK inc Cover PVC Blue 4.8mm Eclipse II V2 | 1009980096 |




- | | | |
|-----|--|------------|
| 12 | SWITCH On/Off Mini/ Eclipse II/ Elite / Admiral II | 1009980053 |
| 13 | FUSEHOLDER PTF35 | 1009020013 |
| 13a | FUSE 1 Amp | 6003990102 |




- | | | |
|----|---|------------|
| 14 | BUTTON Alarm Mute Eclipse II | |
| 15 | CONNECTOR CPC 1/4" BSPT Male Threaded Pump | 1009980085 |
| 16 | CONNECTOR CPC Female Threaded Pump 1/4" BSPT | 1009980086 |
| 17 | SIDE PANEL inc connectors & tubes Prelude / Eclipse / Elite | 1009980088 |




- | | | |
|----|---|------------|
| 18 | ATTENUATOR for Eclipse Pump | 1009980004 |
| 19 | FILTER Ø30mm- PHP1000 Rear Case/PHP398Attenuator | 1009980028 |
| 20 | FEET Back Foam Circular for Attenuator Eclipse Pump | 1009100009 |
| 21 | ATTENUATOR Eclipse Rubber Seal | 1009100006 |




- 22 SEAL Top Hat Rubber Ring Black for Eclipse/Eclipse II Pumps 1009980077
- 23 BELLOW & Internal Housing Compressor Eclipse I/II /Eli 1009980006



- 24 COMPRESSOR Air Dr Assy 3rd Ed (Elite / Eclipse I/II / Elite) 1009980015




- 25 TIMER inc Rotating Valve for Eclipse I/II Pump 1009980005
- 26 FILTER Timer for Eclipse Pump 1009100004
- 27 FILTER CAP Rubber for Eclipse timer 1009100012



28 PCB Main for Eclipse II / Dynabest Pump

1009060001




29 LP Unit for Eclipse Pump

1009100005


30 BAG Low Pressure Sensor

6001020010



31 LOOM Switch LP Eclipse II / Elite Pumps

1009980038




32



33

- | | | |
|----|--|------------|
| 32 | LEAD Mains UK Plug 3rd Ed 2m Straight Grommet w/o Term block (without earth cable) | 1009980007 |
| 33 | LEAD Main UK Plug - Grommet Adl II/Eclip II/Elite /Phase III (with earth cable) | 1009980055 |



34

35

- | | | |
|----|--|------------|
| 34 | CELL Eclipse II / Elite Active Nylon 5" x 85cm Dark blue | 1009980001 |
| 35 | BASESHEET Mattress Eclipse II | 1009060002 |




36



37

- | | | |
|----|--|------------|
| 36 | FEEDPIPE Male for Eclipse II / Phase III | 1009060016 |
| 37 | FEEDPIPE Female for Eclipse II / Phase III | 1009060015 |



38 CPR Pull Tag Assembly 12" Long

1009980058



39 COVER Elite Active System 74" Long Green ▲
COVER Elite Active System 79" Long Green ▲

1004040001
1004040002



40 COVER Eclipse II Active System Jazz (W/O Flaps)
COVER Eclipse II Active System Jazz (With Flaps)
COVER Eclipse II Active System Jazz (W/O Flaps) 80cm Wide

1004060001
1004060002
1004060008