

Park House
Healthcare
By *drive*



2016

PUMP1000 Parts List



L. Parry (Technical)

Park House Healthcare Ltd

28 Jan 2016



- 1 MEMBRANE Front Panel Display only for Eclipse / Pump1000 1009980022
- 2 CASE Front PUMP1000 inc Membrane & Control PCB 1009980018




- 3 CASE Rear Pump1000 inc. Bedhooks, filter & AntiRattle Foam 1009980016
- 4 FOAM Black Strip Anti Rattle Pump 1000 L18cmxW1.5cmxH0.8cm 6003990084
- 5 FUSEHOLDER PTF35 1009020013



- 6 FILTER Ø30mm- PHP1000 Rear Case/PHP398Attenuator 1009980028
- 7 COVER Filter Case Plastic for Pump1000 1009980020



- | | | |
|----|--------------------------------------------------------|------------|
| 8 | BUZZER Alarm 3-24V/3.8KHz 95 dB PHP1000 | 1009980092 |
| 9 | SWITCH On/Off Phase III / Pump1000 | 1009980049 |
| 10 | TIMER Inc Motor & Micro switch Assembly Delta/PUMP1000 | 1009210005 |
| 11 | VALVE Non Return PUMP1000/Phase III 1/4" | |
| 12 | COMPRESSOR Unit (Full) for Pump1000 | 1009980019 |



- | | | |
|----|------------------------------------------------|------------|
| 13 | MEMBRANE Front Panel Display only for Pump1000 | 1009980022 |
|----|------------------------------------------------|------------|



- | | | |
|----|-----------------------------------------------------|------------|
| 14 | BELLOW 230V 50Hz Pumping Set for Phase III/PUMP1000 | 1009980027 |
| 15 | GASKET Bellows 2 holes Black for Phase III/PUMP1000 | 1009010012 |

16 →




16 PCB Main inc EPROM Chip Programmable for Pump1000 1009980025

17 →




17 PCB Control Circuit for Pump1000 1009980024

18 →



18 CABLE Ribbon main PCB to Front Panel Elite Pump1000 1009040010

19 →



← 19

19 FUSE 1 Amp 6003990102
20 FUSE 3 Amp 6003990005




21



22

- 21 LEAD Mains UK Plug 3rd Ed 2m Straight Grommet w/o Term block (without earth cable) 1009980007
- 22 LEAD Main UK Plug - Grommet Adl II/Eclip II/Elite /Phase III (with earth cable) 1009980055



- 23 CONNECTOR CPC 1/4" Male Threaded Pump 1009980085
- 24 CONNECTOR CPC Female Pump 1/4" BSPT 1009980086



25

- 25 CELL Eclipse II / Elite Active Nylon 5" x 85cm Dark blue 1009980001



26



27

- 26 FEEDPIPE Male for Eclipse II / Phase III
- 27 FEEDPIPE Female for Eclipse II / Phase III

1009060016
1009060015



28

- 28 BASESHEET Mattress Eclipse II

1009060002



29

- 29 CPR Pull Tag Assembly 12" Long

1009980058

30



30 COVER Elite Active System 74" Long Green
COVER Elite Active System 79" Long Green

1004040001
1004040002

31



31 FOAM Elite Active System Insert Visco 74"x35"x2" Pink Elite
FOAM Elite Active System Insert Visco 79"x35"x2" Pink Elite

1009040001
1009040002

32



32 COVER Eclipse II Active System Jazz (W/O Flaps)

1004060001