



- | | | |
|--|---|--|
| <ol style="list-style-type: none"> 1. SPEED WIGWAG
CODE: HW-72011093 2. THROTTLE ASSEMBLY
CODE: HW-77126120 3. HANDLE GRIP
CODE: HW-78042039 4. BUZZER
CODE: HW-77126025 5. HANDLE BAR
CODE: HW-32052234 6. PLUG OF HANDLE BAR
CODE: HW-72012006 7. PLUG OF STEM
CODE: HW-72012005 8. LOCKING KNOB
CODE: HW-78132053 9. HANDLE TUBE WELDING
CODE: HW-32052244 10. HEADLIGHT
CODE: HW-77196075 11. HEADLIGHT COVER
CODE: HW-77192102 12. HORN BUTTON
WITH WIRING
CODE: HW-77126009 13. LITHIUM BATTERY
CHARGER
CODE: HW-77086041 14. CHARGING SOCKET
WIRING
CODE: HW-77146129 15. UPPER BODY COVER
CODE: HW-71272774, HW-71272791, HW-71272793 16. GYROSCOPE WIRING
CODE: HW-77122010 17. FRONT LOWER COVER
CODE: HW-71272755 18. ENGINE LIGHT EXTENSION
CORD SET
CODE: HW-77206132 19. MAIN WIRING
CODE: HW-77156094 20. SLIDING TUBE
CODE: HW-31161093 21. IGNITION KEY SOCKET
WIRING,
W/2 KEYS
CODE: HW-77146107 22. HEADLIGHT SWITCH
WIRING
CODE: HW-77146040 23. SPEED ADJUSTING KNOB
CODE: HW-77126024 24. HANDLE BAR UPPER
COVER
CODE: HW-72042415,
HW-72042426,
HW-72042427 25. LED AMMETER LI BATTERY
CODE: HW-77206139 26. SPEED POTENTIOMETER
WIRING
CODE: HW-77146017 27. UPPER CONTROL PC BOARD
CODE: HW-77126075 | <ol style="list-style-type: none"> 28. BATTERY BOX COVER
CODE: HW-77172577, HW-77172590, HW-77172591 29. BODY PROTECTION COVER
CODE: HW-71273002 30. OVAL STICKER
CODE: HW-70010708 31. BATTERY WIRE
CODE: HW-77186162 32. LITHIUM BATTERY
CODE: HW-77172581 33. ELECTRIC FOLDING
CONTROL
CODE: HW-77107200 34. ACTUATOR
CODE: HW-73066044 35. BOTTOM PIPE WELDING
CODE: 3HW-1162223 36. ELECTRIC FOLDING
REMOTE CONTROL
CODE: HW-77107212, HW-77107216, HW-77107217 37. ELECTRIC FOLDING
CONTROL
LINE GROUP
CODE: HW-77156095 38. METAL BUTTON
WITH GREEN LIGHT
CODE: HW-77206135 39. AFTER THE
DECORATIVE COVER
CODE: HW-71272761 40. METAL BUTTON
WITH RED LIGHT
CODE: HW-77206134 41. TAILLIGHT
CODE: HW-77196076 42. R PU REAR WHEEL
CODE: HW-74055063 43. TAILLIGHT & GYROSCOPE
EXTENSION WIRING
CODE: HW-77206133 44. CONTROLLER
CODE: HW-77107207 45. ANTI-TIPPING WHEEL
CODE: HW-74232014 46. TRANSAXLE
CODE: HW-77017567 47. FRAME WELDING
CODE: HW-31162183 48. L PU REAR WHEEL
CODE: HW-74055043 49. WHEEL CAP
CODE: HW-74172082 50. MAT
CODE: HW-71284036 51. FRONT UPPER COVER
CODE: HW-71272773,
HW-71272792,
HW-71272794 52. FRONT FRAME WELDING
CODE: HW-31152203 53. STEERING AXLE WELDING
CODE: HW-31151183 54. FRONT WHEEL COVER
CODE: HW-74212005 | <ol style="list-style-type: none"> 55. PU FRONT WHEEL
CODE: HW-74055042 56. PUSH TUBE WELDING
CODE: HW-31152213 57. ROTATING INNER TUBET
CODE: HW-71401128 58. ANTI-TIPPING WHEEL
CODE: HW-74232015 59. LOWER BODY COVER
CODE: HW-71272759 60. FRONT POST SLIDE WELDING
CODE: HW-33042434 61. OUTER CROSS
BRACE WELDING
CODE: HW-33042414 62. INNER CROSS
BRACE WELDING
CODE: HW-33042404 63. POSITIONING ROD
CODE: HW-71401169 64. SEAT ASS'Y
CODE: HW-73054098 65. SEAT INNER FRAME WELDING
CODE: HW-33042463 66. SEAT OUTER
FRAME WELDING
CODE: HW-33042463 67. REAR ROTATING INNER TUBE
CODE: HW-71401129 68. ROTATING SPACER
CODE: HW-73012037 69. HANDLE BAR LOWER COVER
CODE: HW-72042406 70. POWER CABLE FOR CHARGER
CODE: HW-77096015 71. PLUG
CODE: HW-78022061 |
|--|---|--|