

## **PRODUCT INSTRUCTION SHEET**

### **PRODUCT CODE**

Foxton Shower Chair Fixed Ht with Handle (140)  
Foxton Shower Chair Fixed Ht without Handle (140NHANDLE)

### **DESCRIPTION**

Nylon coated all steel welded frame with 'H' section base for added strength. Moulded polypropylene seat with drainage holes. Available with or without handle.



### **TECHNICAL SPECIFICATIONS**

Internal width between arms	18" (140) & 19" (140A)
External width:	20" (140)
Seat Width x Depth:	17" x 15"
Floor to seat height	20"
Overall product height:	32.5"
Product Weight:	6.2kg
Weight Capacity:	165kg (26st)

### **USAGE REQUIREMENTS**

1. OTs should do a risk assessment on any appliance at time of issue, ensuring the products surface stability.
2. Do not exceed weight limit of 165kg (26 Stone) and ensure that patients do not 'swing-on' the frames either in a forward or backward motion.
3. Clean frame with cream cleaner, seat with soap and water.
4. This product is intended for use in a wet room / shower. It is not intended for use outside and should not be immersed in water for long periods of time. Ensure the product is dried thoroughly after use.

### **WARRANTY**

Your Drive DeVilbiss branded product is warranted to be free of defects in materials and workmanship for one year from date of purchase. This device was built to exacting standards and carefully inspected prior to shipment. This warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the customer. In the event of a defect covered by this warranty, we will at our option supply parts or replace the device. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, castors and grips, which are subject to normal wear and need periodic replacement. If you have any queries or require full warranty conditions, please contact your Drive DeVilbiss stockist. This does not affect your statutory rights.