

Crystallography materials and methods for PTP1B-Compound **8e** and **22i**.

Co-crystallization of PTP1B with inhibitors

A 6-10 mg/mL preparation of PTP1B in 10 mM Tris pH 7.5, 25 mM NaCl, 0.2 mM EDTA and 3 mM DTT, was used for crystallization. Crystals were grown by the sitting, as well as the hanging drop, vapor diffusion methods. A 1:10 (PTP1B:inhibitor) molar ratio mixture was prepared at least one hour prior to crystallization. Two μL PTP1B-inhibitor solution was mixed with 2 μL reservoir solution consisting of 0.1 M HEPES buffer pH 7.5, 0.3-0.4 M Na-acetate or Mg-acetate, 12-16% (w/v) polyethylene glycol 8000 and/or 4% (v/v) glycerol. The reservoir volume was 1 mL. Crystals grew to the size of 0.3-0.6x0.1-0.3x0.1-0.3 mm within 2-3 days.

Data collection.

All diffraction data collections were performed at 100 K. The following cryo conditions were used: to the hanging or sitting drop 3 μL of 50% (v/v) glycerol (containing 0.5 mmol inhibitor) were added. The crystal was removed from the drop after 5-30 min and transferred to 50% glycerol (containing 0.5 mmol inhibitor) and flash frozen. Data were collected using a Mar345 image plate detector at the MAX-lab synchrotron facilities at Lund University. A 1° oscillation per image was used for 60 images. The space group was determined to be P3121 for the crystals obtained with compound **8e** and P43212 for all the crystals obtained with compound **22i**. Data processing was performed using Denzo, Scalepack and the CCP4 program suite.^{1,2}

Refinements.

As P3121 contains a polar axis and, thus, possesses more than one indexing possibility, a molecular replacement solution using Amore³ was determined prior to the refinements. Further, a molecular replacement solution was also calculated for the P43212 space group. A high resolution PTP1B structure was used as a starting model (PDB file: 1PTV),⁴ with ligand and water molecules omitted from the structure. All refinements were performed with CNX (Asselrys Inc.). Interchanging cycles of model building using Quanta X-build (Asselrys Inc.) and refinement were performed. The $2F_o - F_c$ maps were inspected by the use of X-ligand (Asselrys Inc.) at a 1 sigma level for densities that could correspond to the structures of the inhibitors. In both cases, a well-suited inhibitor electron density was identified in the active site pocket (see below). No other densities were identified to fit the inhibitors. Water molecules were inserted using the X-solvate program (Asselrys Inc.) based on 1.5σ $2F_o - F_c$ electron density maps. For further details see Table 1.

References

1. Otwinowski, Z. and Minor, W. Processing of X-ray diffraction data collected in oscillation mode. In *Methods Enzymol*, ed.; Carter, C. W. and Sweet, R. M., Eds. Academic Press, New York, 1997; vol. 276, pp 307- 326.
2. Collaborative Computational Project, N. 4. The CCP4 Suite: Programs for protein crystallography. *Acta Crystallogr., Sect. D* 1994, 50, 760-763.
3. Navaza, J. Amore - an automated package for molecular replacement. *Acta Crystallogr., Sect. A* 1994, 50, 157-163.

4. Jia, Z. C., Barford, D., Flint, A. J., and Tonks, N. K. Structural basis for phosphotyrosine peptide recognition by protein tyrosine phosphatases 1B. *Science* **1995**, *268*, 1754-1758.

Table 1. Statistics of X-ray data and structure refinements.

Table 1.

| Protein | PTP1B | PTP1B |
|--|---------------------|---------------------|
| Ligand | Compound 8e | Compound 22i |
| Space group | P3121 | P43212 |
| Unit cell parameters | a=b=88.2 c=104.1 Å | a=b=97.74 c=92.347Å |
| Completeness (all data) | 98.9% (30-1.89 Å) | 94.3 % (30-2.0 Å) |
| Completeness (high res.) | 99.8% (1.91-1.89 Å) | 98.8% (2.03-2.00 Å) |
| Multiplicity (all data) | 3.5 (30-1.89 Å) | 4.4 (30-2.0 Å) |
| Rmerge (all data) | 8.2% (30-1.89 Å) | 4.3% (30-2.0 Å) |
| Rmerge (high res.) | 49% (1.91-1.89 Å) | 44.7% (2.03-2.00 Å) |
| <I/σ(I)> (all data) | 14.3 (30-1.89 Å) | 20.7 (30-2.0 Å) |
| <I/σ(I)> (high res.) | 2.7 (1.91-1.89 Å) | 2.2 (2.03-2.00 Å) |
| Unique reflections | 37618 | 29166 |
| Atoms in structure | 2764 | 2683 |
| R-factor ¹ | 20.5% | 22.7% |
| R-free ² | 21.9% | 28.5% |
| R.m.s. deviations from idealised geometry's | | |
| Bond lengths (Å) | 0.006 | 0.006 |
| Bond angles (°) | 1.2 | 1.3 |

¹Crystallographic R-factor – calculated using all data. ²Rfree = $\frac{\sum_{(hkl) \in T} |F_o| - |F_c|}{\sum_{(hkl) \in T} |F_o|}$

where T is a test set containing a random 5% of the observations omitted from the refinement process.