

Introduction

- Collapse is referred as the loss of the pore structure of freeze-dried cake whereas “**Micro-collapse**” is an intact cake with increasing pore in dried layer and can promote drying.



- Dry product resistance (R_p)** is one of critical variables having a great impact on primary drying process.
- R_p is defined as the resistance to mass flow of water vapour from the product through the pores structure in the dried layer and can be expressed by:

$$\frac{dm}{dt} = A_p \frac{\Delta P}{R_p}$$

where $\frac{dm}{dt}$ is the drying rate (g/h/vial), A_p is the internal cross-sectional area of the vial, ΔP is the pressure difference between pressure of ice at sublimation interface and chamber pressure and R_p is the area normalized resistance of the dried product.

- An increase in product temperature during primary drying stage above the **Collapse Temperature (T_c)** may cause the collapse of a freeze-dried cake with the possible rejection of the entire production batch
- Collapse temperature could be determined by
 - A freeze-drying microscope (FDM)
 - An optical coherence tomography based freeze drying microscopy (OCT-FDM), Mujat (2012)
- However, these current techniques have some limitations
 - FDM → Off-line measurement, Not always provide the information of micro-collapse, Varshney (2015)
 - OCT-FDM → Single vial technique
- Through vial impedance spectroscopy (TVIS)**, a novel non-invasive techniques has been shown previously to be sensitive to the collapse event itself, through dramatic changes in the electrical capacitance of a solution filled in freeze-drying vial, Smith (2014).

Aim

- To evaluate the application of TVIS system for the prediction of micro-collapse during a freeze-drying cycle.

Materials and Methods

- The electrical impedance of a 5%w/v lactose solution contained within modified glass freeze-drying vial (TVIS vial, Fig.1A) was measured over the frequency range of 10 Hz to 1 MHz by using TVIS system during a freeze-drying process.
- A full load of vials with TVIS vial at the center (Fig.1B) was then placed on a single shelf of a Virtis Advantage Plus benchtop Freeze-dryer.
- A freeze drying protocol with an annealing step is performed to dry the solution. A critical feature of the drying cycle is the inclusion of a ramp in the shelf temperature during primary drying that will force the product through its collapse event.
- Scanning electron microscopy (SEM) images of the freeze-dried cake were acquired at a 500x magnification.




Fig. 1. (A) TVIS vial (left) and a neighboring vial with thermocouple (right), (B) The cluster of vials with TVIS vial at center

Results and Discussions

- The correlation between $\text{Log} F_{PEAK}$ from TVIS vial and the thermocouple temperature in a neighboring vial of the re-heating part of annealing step (Fig.2A) provides a predictive product temperature at primary drying process defined as **T-FPEAK** (Fig.2B).
- The decrease in C_{PEAK}'' parameter corresponding to the amount of ice bounded within electrode region can be used to estimate drying rate. However, this parameter also depends on temperature. As the temperature is increased, C_{PEAK}'' value increases (Fig.2C). Therefore, the temperature compensation for this parameter defined as normalized C_{PEAK}'' (\hat{C}_{PEAK}'') is required (Fig.2D) by using the **temperature correction factor (ϕ)** which is calculated from the re-heating phase of the annealing stage and C_{PEAK}'' at starting temperature of primary drying as the reference value.
- C_{PEAK}'' can be normalized by the following equation:

$$\hat{C}_{PEAK}'' = \frac{C_{PEAK}''(T_t)}{\phi}$$

where $C_{PEAK}''(T_t)$ is C_{PEAK}'' at time (t) and temperature (T) during primary drying and ϕ is the temperature correction factor from C_{PEAK}'' calibration of re-heating step




Fig. 2. The temperature and TVIS parameters profile of 5%w/v lactose solution during the primary drying stage. (A) Temperature calibration from re-heating step, (B) a predictive temperature during primary drying, (C) C_{PEAK}'' calibration from re-heating step and (D) Temperature-compensated C_{PEAK}'' (\hat{C}_{PEAK}'')

- At 4.8 hour into primary drying there is a significant increase in the rate of change of \hat{C}_{PEAK}'' which corresponds to an dramatic increase in drying rate as shown in Fig.3B. This suggests there is a microscopic change in cake structure, due to micro-collapse, which results in an increase the pore size distribution in the freeze-dried matrix thereby decreasing the product resistance (Fig.3C) and consequently improving vapour flux. The predicted temperature at this point in time is equal to the collapse temperature of -32°C (Fig.3A).
- This suggestion is confirmed by cake morphology images by SEM as shown in Fig.3E. At dried layer thickness of 0.27 cm corresponding to 4.8 hour of primary drying Fig.3D, a micro-collapse layer has developed which can be demonstrated by SEM as larger pores in middle layer.




Fig. 3. Results of 5%w/v lactose solution during the primary drying stage (A) The temperature (T-FPEAK) profile, (B) \hat{C}_{PEAK}'' parameters and drying rate over drying period, (C) Product resistance (R_p) as a function of drying time, (D) Product resistance (R_p) as a function of dry layer thickness and (E) SEM of top and middle layer of lactose cake at the end of the cycle.

Conclusions

- A significant decrease in C_{PEAK}'' at the point of micro-collapse (as confirmed by SEM) highlights the potential for using TVIS for monitoring microscopic changes within cake during primary drying step.
- This study demonstrates a prospective use of TVIS as a process control tool that would allow the cycle to be driven at the highest achievable temperature whilst avoiding collapse.

References

- Mujat, M., Greco, K., Galbally-Kinney, K.L., Hammer, D.X., Ferguson, R.D., Iftimia, N., Mulhall, P., Sharma, P., Pikal, M.J., Kessler, W.J., 2012. Optical coherence tomography-based freeze-drying microscopy. *Biomed. Opt. Express*, 3, 55-63.
- Smith, G., Arshad, M.S., Nazari, K., Polygalov, E., Ermolina, I., Taylor, J., Page, T., 2014. Through-Vial Impedance Spectroscopy: A New In-Line Process Analytical Technology for Freeze Drying. *Pharm Technol*, 38.
- Varshney, D., Singh, M., 2015. *Lyophilized Biologics and Vaccines: Modality-Based Approaches*, Springer, New York.

Acknowledgements

GEA Pharma Systems and AstraZeneca (2008-12) and co-funded by Innovate UK (formerly the Technology Strategy Board).

