



SPECIALITIES


**Intermaxillary Fixation
Micromotor System
Alveolar Bone Graft
Kretschmer Reference Pin
Further Products**



**Intermaxilläre Fixation
Intermaxillary Fixation
Fijación intermaxilar
Fixation intermaxillaire
Fissaggio intermascellare**


“POWER GRIP”

Schraubenkopf mit kreuzförmigem Schlitz
Screw head with cross slot
Tornillo de cabeza ranurada en cruz
Tête de vis avec fente cruciforme
Testa della vite a croce



Verpackungseinheit 5 Schrauben im Magazin (nicht steril)
Packing unit 5 screws in storage tray (non-sterile)
Paquete de 5 tornillos en cartucho (no estéril)
Emballage unitaire de 5 vis en cartouche (non stériles)
Confezione da 5 viti nella cartuccia (non sterile)

Festziehen des Drahtes über die Taille
Tightening of the wire over the waist
Fijación del alambre sobre el cuello
Serrage du fil au-dessus du cou
Fissaggio del filo d'acciaio attraverso il collo



Durchführen des Drahtes durch die Bohrungen
Pulling of the wire through the drillings
Pasar el alambre por los orificios
Filetage du fil par les percements
Passare il filo d'acciaio tramite le perforazioni



Bei harter Kortikalis oder Fixation kleiner Knochenfragmente wird ein Vorbohren empfohlen.
In case of hard cortical bone substance or fixation of small bone fragments we recommend pre-drilling.
En caso cortical dura o de la fijación de fragmentos pequeños recomendamos pretaladras.
Pour le cas où la corticale serait trop dure ou pour la fixation de petits fragments osseux, nous recommandons un préforage.
In presenza di una corticale molto dura oppure per la fissazione di piccoli frammenti consigliamo di perforare prima dell'inserimento della vite.

INTERMAXILLARY FIXATION

Titan-Schrauben
Titanium screws
Tornillos de titanio
Vis en titane
Viti in titanio



 P			
	a	Self-drilling Ø 2,0 mm	Self-tapping Ø 2,0 mm
	5,0 mm	68.94.29B	
	8,0 mm	68.94.28B	68.94.30B
	10,0 mm	68.94.26B	68.94.31B
	12,0 mm	68.94.27B	68.94.32B
	14,0 mm		68.94.33B
	16,0 mm		68.94.34B

Schraubendrehergriff
Screwdriver handle
Mango para atornillador
Manche de tournevis
Manico per cacciavite



		10 cm - 4"	75.23.98A	2.0 2.4
--	--	------------	------------------	-----------------------


Schraubendreherklingen
Screwdriver blades
Hojas de atornillador
Lames de tournevis
Lame per cacciavite



	"POWER GRIP"	9 cm - 3 1/2"	75.23.00	2.0
--	---------------------	---------------	-----------------	------------



mit Greifvorrichtung
with gripping device
con garra
avec préhenseur
con dispositivo da presa



	"POWER GRIP"	9 cm - 3 1/2"	75.23.40	2.0
--	---------------------	---------------	-----------------	------------



Selbstbohrend
Self-drilling
Autoperforante
Autoperforante
Autoperforante




Selbstschneidend
Self-tapping
Autoroscante
Autotaraudeuse
Autofilettante



Intermaxilläre Fixation
Intermaxillary Fixation
Fijación intermaxilar
Fixation intermaxillaire
Fissaggio intermascellare

Spiralbohrer
 Twist drills
 Brocas espirales
 Forets hélicoïdaux
 Perforatori elicoidali




	5 - 10 mm	a = 11 mm	b = 44 mm	68.76.73	☐ ☐ 2.0 ☐
	12 - 16 mm	a = 18 mm	b = 50 mm	68.76.88	☐ ☐ 2.0 ☐

Flachzange
 Flat nose pliers
 Alicata plano
 Pince plate
 Tenaglia piatta



			14 cm - 5½"	14.46.14	☐ ☐ 2.0 ☐
--	--	--	-------------	-----------------	---

Drahtschneidezange
 Wire cutting pliers
 Alicata para cortar alambre
 Pince coupe-fils
 Pinza per tagliare fili



			12 cm - 4¾"	14.30.11	☐ ☐ 2.0 ☐
--	--	--	-------------	-----------------	---

Drahtspannzange
 Wire twisting forceps
 Pinza tensora de alambre
 Pince serre-fils
 Pinza per tendere fili



	Obwegeser		13 cm - 5½"	68.04.73	☐ ☐ 2.0 ☐
--	------------------	--	-------------	-----------------	---



INTERMAXILLARY FIXATION

Drahtinstrumente
Wire instruments
Instrumentos para alambre
Instruments pour fils
Strumenti per fili



Luniatschek	 a = 2 mm 	18 cm - 7 1/8"	330.01.86	☐ ☐ 2.0 ☐
Luniatschek	 a = 3 mm 	18 cm - 7 1/8"	330.01.80	☐ ☐ 2.0 ☐
Luniatschek	 a = 3 mm 	18 cm - 7 1/8"	330.01.85	☐ ☐ 2.0 ☐

Drähte
Wires
Alambres
Fils
Fili





	Ø = 0,3 mm (29 G)	10 m	76.01.03	☐ ☐ 2.0 ☐
	Ø = 0,4 mm (27 G)	10 m	76.01.04	☐ ☐ 2.0 ☐
	Ø = 0,5 mm (25 G)	10 m	76.01.05	☐ ☐ 2.0 ☐

Kassette (ohne Implantate, Schrauben und Werkzeuge)
Tray (without implants, screws and tools)
Caja (sin implantes, tornillos y sus respectivos instrumentos)
Cassette (sans implants, vis et outils)
Cassetta (senza impianti, viti e strumentario)



IMF Intermaxillary Fixation	23 x 11,5 x 4 cm	69.02.26	☐ ☐ 2.0 ☐
--	------------------	-----------------	---





Intermaxilläre Fixation
Intermaxillary Fixation
Fijación intermaxilar
Fixation intermaxillaire
Fissaggio intermascellare


IMF SET



Satz komplett
Set complete
Juego completo
Jeu complet
Set completo

IMF Intermaxillary Fixation	23 × 11,5 × 4 cm	68.77.92	☐ ☐ 2.0 ☐
--	------------------	-----------------	---

Bestehend aus:
Consisting of:
Compuesto por:
Composé de:
Composto da:



01	68.76.73	Twist drill	1x
02	68.76.88	Twist drill	1x
03	68.94.26B	Screw for IMF, self-drilling, 10 mm, P/5	1x
04	68.94.27B	Screw for IMF, self-drilling, 12 mm, P/5	1x
05	68.94.28B	Screw for IMF, self-drilling, 8 mm, P/5	1x
06	68.94.29B	Screw for IMF, self-drilling, 5 mm, P/5	1x
07	68.94.30B	Screw for IMF, self-tapping, 8 mm, P/5	1x
08	68.94.31B	Screw for IMF, self-tapping, 10 mm, P/5	1x
09	68.94.32B	Screw for IMF, self-tapping, 12 mm, P/5	1x
10	68.94.33B	Screw for IMF, self-tapping, 14 mm, P/5	1x
11	68.94.34B	Screw for IMF, self-tapping, 16 mm, P/5	1x
12	75.23.00	Screwdriver blade	1x
13	75.23.98A	Screwdriver handle	1x
14	69.02.26	Tray (without implants, screws and tools)	1x

P_X Verpackungseinheit "X" Stück
Packing unit "X" pieces
Paquete de "X" unidades
Emballage de "X" pièces
Confezione da "X" pezzi

