

MagniLink PCViewer™

– innovative software for PC

For increased participation -
in school, at work
and at home!



MagniLink PCViewer™

MagniLink PCViewer™ is an innovative software for Windows based PC computers that is used to display the image from a MagniLink reading camera system. The software controls the reading camera's features, such as magnification, artificial colors and brightness, and provides the user with the ability to make individual adjustments for the best possible experience.

Functionality of MagniLink PCViewer™:

- Different color and brightness settings, e.g., natural and artificial colors.
- Different levels of magnification.
- Lockable autofocus.
- Take and view pictures.
- Record and view videos.
- Mirror, rotate and freeze pictures.
- Automatic updates (requires Internet).
- Various guide lines to facilitate reading.
- Now possible to change and adapt the product even more in EUC tab (End-UserConfiguration).
- Different ways to control functions. Either via menus, keyboard shortcuts or buttons on the reading camera's control panel.
- Configuration of own shortcuts.
- Camera image positioned at top/front of all open applications.
- Select several active views and switch-between them.
- Compatible with the ZoomText and SuperNova programs.

Options:

The TTS option makes it possible to OCR process text and have it read aloud using integrated speech synthesis. The OCR processed text can be displayed in several different picture view modes synchronized with the speech. It is possible to OCR process text from a picture taken with the connected MagniLink camera, and it is also possible to import other documents into the program to have them OCR processed. It is possible to save OCR result



as a sound file (.wav) and you can enlarge and navigate in the OCR-picture. It is also possible to paste text into the program.

Recommended system requirements:

- Microsoft Windows 7.
- 1 integrated USB 3.0-port
- i5 or i7 processor, 4 GB RAM, external graphics card with at least 512 MB built-in memory.

General:

MagniLink PCViewer™ is designed for PC computers. There is similar software for Chromebook and Mac computers.



LVI Low Vision International
Verkstadsgatan 5,
352 46 Växjö, SWEDEN
Tel: +46 (0)470 72 77 00
Fax: +46 (0)470 72 77 25
E-mail: info@lvi.se
WWW.LVI.SE