

MagniLink PCViewer™

– innovative software for PC

For increased participation -
in school, at work
and at home!



MagniLink PCViewer™

MagniLink PCViewer™ is an innovative software for Windows based PCs that is used to display the image from a MagniLink video magnifier. The software controls the reading camera's features, such as magnification, artificial colors and brightness, and provides the user with the facility to make individual adjustments for the best possible experience.

Functionality of MagniLink PCViewer™:

- Different color and brightness settings, e.g., natural and artificial colors.
- Different levels of magnification.
- Lockable autofocus.
- Take and view pictures.
- Record and view videos.
- Mirror, rotate and freeze pictures.
- Automatic updates (Internet).
- Various guide lines to facilitate reading.
- Now possible to change and adapt the product even more in EUC tab (End User Configuration).
- Different ways to control functions. Either via menus, keyboard shortcuts or buttons on the reading camera's control panel.
- Configuration of own shortcuts.
- Camera image positioned at top/front of all open applications.
- Select and switch between several views.
- Compatible with the ZoomText and SuperNova programs.

Options:

The TTS option makes it possible to OCR process text and have it read aloud using integrated speech synthesis. The OCR processed text can be displayed in several different picture view modes synchronized with the speech. It's possible to OCR process text from a picture taken with the connected MagniLink camera, and it is also possible to import other documents into the application to have them OCR processed. You can enlarge and navigate in the OCR-picture and the OCR result can be saved as an audio file (.wav). It's also possible to paste text into the application.



If you purchase the optional TTS, you can also use the software standalone ie without a connected video magnifier. Then you can get already saved text or image documents OCR processed and the text read aloud using speech synthesis. When purchasing the video magnifier MagniLink S Premium 2 or MagniLink Pro, the software PCViewer is included and if you also choose the option TTS the standalone feature is also included. You can also buy TTS and standalone afterwards.

Recommended system requirements:

- Microsoft Windows 7.
- 1 integrated USB 3.0-port
- i5 or i7 processor, 4 GB RAM, external graphics card with at least 512 MB built-in memory.



LVI America, Inc
150 N. Michigan Ave., Suite 1950
Chicago, Illinois 60601, USA
Phone: (888) 781-7811
E-mail: info@lviamerica.com
WWW.LVIAMERICA.COM