

Let us blow you away with
a brand new software for Mac!

MagniLink MacViewer

We are pleased to announce MagniLink MacViewer – a software that is used to display image from a MagniLink video magnifier on a Mac. The software controls features such as magnification, color and brightness, and provides the user with the facility to make individual adjustments for the best possible experience. Enjoy!



This is what you get:

- Different color and brightness settings, e.g., natural and artificial colors.
- Different levels of magnification.
- Lockable autofocus.
- Take and view pictures.
- Record and view videos.
- Mirror, rotate and freeze pictures.
- Various guide lines to facilitate reading.
- Possible to change and adapt the product even more in EUC tab (EndUserConfiguration).
- Different ways to control functions. Either via menus, keyboard shortcuts or buttons on the reading camera's control panel.
- Configuration of own shortcuts.
- Camera image positioned at top/front of all open applications.
- Select several active views and switch between them.
- And several other smart fetures!

Options:

The TTS option makes it possible to OCR process text and have it read aloud. The OCR processed text can be displayed in several different picture view modes synchronized with the speech. It is possible to OCR process text from a picture taken with the connected MagniLink camera, and it is also possible to import other documents into the program to have them OCR processed.



Interested?
– just contact
us and we tell
you more!



LVI America, Inc
302 Saunders Road, Suite 200, Riverwoods, IL 60015, USA
Phone: (888) 781-7811 E-mail: info@lviamerica.com

LVI assumes no responsibility for any inaccuracies that may be contained in this sheet.