

MagniLink Voice 2



Give your text a voice!

MagniLink Voice 2 is a reading machine with focus on simplicity, high performance and contemporary design. It's the user-friendly reading machine that exceeds expectations!

Easy to use

MagniLink Voice 2 is extremely easy to use – just push the orange button on the control panel and convert the text into speech!

Design and function in harmony

Premium material, details and finish ensure that the unit is as stable as it is durable. Robust, yet extremely compact, it is designed for traveling and moving around. And when placed at home or office, it looks great in any environment.

Adapt with options

If you choose a unit with monitor output, you can present the scanned text on a monitor, thus adding a visual display to the text being read aloud. This facilitates

document orientation and improves reading comprehension. If you choose a reading machine with a battery you can listen to your text wherever you are!

Other features

- Easy set up – just move the camera arm up and start.
- Adaptable to your specific needs, by adjusting reading speed, volume, audio profiles, voices, light and more.
- Choose the standard settings mode or the advanced settings mode.
- Automatic language detection.
- Extremely good depth of field – optimized for browsing large books.
- Documents on USB can also be read aloud.
- Voice feedback for all functions.
- Reads a full A4/US Letter page.





Technical specifications

Model	MagniLink Voice 2
Dimensions	
Folded:	23x6.5x32,5 cm (9x2.5x12.8 in)
Unfolded:	23x31x32 cm (9x12.2x12.6 in)
Weight	2.4 kg (5.3 lbs)
Battery (optional)	
Operating time:	5 h
Charging time:	2.5 h



LVI Low Vision International is recognized as one of the world's leading manufacturers of low vision aids with the mission to make daily life easier for people with visual impairments.

We are certified to ISO 14001 & 9001



Intertek

MagniLink Voice 2 20200318



LVI Low Vision International
Verkstadsgratan 5
352 46 Vxj, SWEDEN
Tel: +46 (0)470 72 77 00
Fax: +46 (0)470 72 77 25
E-mail: info@lvi.se
WWW.LVI.SE