

MagniLink Vision Full HD

Ergonomic and user-friendly

MagniLink Vision is a video magnifier that offers users with visual impairments all the functions they need at an attractive price. It's a modern, flexible video magnifier with an ergonomic and user-friendly design.



It's your choice!
– choose a unit with or without Text To Speech function.

Plug & play

Connect MagniLink Vision to a power socket, press the on/off button and you're ready to start using your video magnifier!

Easy for all

MagniLink Vision is very easy to use, which makes it very popular among elderly users as well.

The control buttons are tactile and ergonomically placed directly on the integrated LVI monitor.

The high-quality monitor is low vision-adapted with high brightness as well as good contrast and color reproduction. The surface of the monitor is matt to avoid reflections.

The high-quality reading table has friction brakes in the X and Y axis. The brakes can be used to provide higher resistance when the reading table is moved sideways and/or forwards/backwards. They can also be used to lock the reading table during transportation.

MagniLink Vision Full HD TTS

The Text To Speech option offers the best of both worlds – a high-quality reading camera and reading machine in one single product, resulting in increased reading speed, reading durability and a more relaxed reading experience.

Please note that TTS must be delivered from the start. It's not possible to order TTS afterwards.

Useful features with TTS

- Read a full US Letter page
- Possible to open/save text on a USB
- Reads misaligned documents
- Automatic language detection
- Adjustable size, color and font
- Listen to the text through built-in speakers or plug in headphones

Technical specifications:

Model	MagniLink Vision MagniLink Vision TTS
Screen size	23"
Magnification	2.0-70x
Resolution	1920x1080
Dimensions in use	22x18.5x22.4 in (WxDxH)
Weight	29.7 lbs

This is included:
Reading camera with a reading table and monitor, polishing cloth, user guide and power adapter. If you have ordered a unit with TTS, you also get a hand control.

With a TTS-unit

– just push the orange button on the hand control and convert the text into speech!



LVI America, Inc
302 Saunders Road, Suite 200, Riverwoods, IL 60015, USA
Phone: (888) 781-7811 E-mail: info@lviamerica.com

LVI assumes no responsibility for any inaccuracies that may be contained in this sheet.