
MagniLink Software



Fits all your needs!

The software is used to display the image from your MagniLink video magnifier, and you choose software based on your computer – PC, Mac or Chromebook. From the software you control the video magnifiers' features, such as magnification, artificial colors and brightness. The software provides you with the ability to make individual adjustments for the best possible experience!

Functionality

Although MagniLink PCViewer, MagniLink MacViewer and MagniLink ChromeViewer are not completely identical, most of the features are the same. For example:

- Different color and brightness settings.
- Different levels of magnification.
- Lockable autofocus.
- Take and view pictures.
- Record and view videos.
- Mirror, rotate and freeze pictures.
- Various guiding lines.
- Different ways to control functions; via menus, keyboard shortcuts or buttons on the reading camera's control panel.
- Configuration of own shortcuts.
- Camera image positioned at top/front of all open applications.
- Select and switch between several views.
- Support for motorized cameras in Magni-Link AIR family.
- Support for multiple video streams and split image between them, e.g. reversal of multiple MagniLink cameras.
- Note tool with the ability to draw, pan and enlarge the image (not MacViewer).

Optional TTS

The TTS (Text To Speech) option for PC, Mac and Chromebook makes it possible to OCR process text and have it read aloud. The OCR processed text can be displayed in several different picture view modes synchronized with the speech. It's possible to OCR process text from a picture taken with the connected MagniLink camera, and it is possible to enlarge and navigate in the OCR-picture. For PC and Mac it's also possible to import other documents into the application or paste texts straight into the application to have them read aloud.

MagniLink PCViewer Standalone

MagniLink PCViewer Standalone is a software that can be used on a Windows based PC without being connected to a MagniLink video magnifier. With the software, you can OCR-process documents and get the text read aloud by the built-in speech synthesis. You can also display the OCR processed text with or without speech in different viewing modes. When purchasing the video magnifier MagniLink S Premium 2 or MagniLink PRO, the "regular" software PCViewer is included and if you choose the option TTS, the standalone feature is also included. MagniLink PCViewer Standalone can also be purchased separately.

Compatibility

MagniLink PCViewer is compatible with other magnification and screen reading software like ZoomText and SuperNova.



Recommended system requirements

MagniLink PCViewer


- Microsoft Windows 7
- One integrated USB 3.0-port
- i5 or i7 processor, 4 GB RAM, external graphics card with at least 512 MB built-in memory

MagniLink MacViewer

- Mac Os 10.14 or later
- One integrated USB 3.0-port
- i5 or i7 processor, 4 GB RAM, external graphics card with at least 512 MB built-in memory

MagniLink ChromeViewer

- One integrated USB 3.0-port
- i5 processor
- 4 GB RAM



MagniLink Software

LVI Low Vision International is recognized as one of the world's leading manufacturers of low vision aids with the mission to make daily life easier for people with visual impairments.

We are certified to ISO 14001 & 9001





LVI Low Vision International Verkstadsgatan 5 352 46 Växjö, SWEDEN Tel: +46 (0)470 72 77 00 Fax: +46 (0)470 72 77 25

E-mail: info@lvi.se

WWW.LVI.SE