



**Comments:** Good correlation with reference instrument which improves when concentrations are above 500ppb due to the poor performance of reference instruments at low concentrations. In yellow color can be appreciated the malfunction of the reference instrument which reduce its availability.

Daily variations are very similar as well as the average concentration.

Version: KUNAKAIR V2.







R <sup>2</sup>	ACCY	90% CONF.	DOQ	AVG. GC.	AVG. GC.	AV
	(PPB)	(PPB)	(PPB)	REF (PPB)	DUT (PPB)	REF/DUT
0,69	86,79	177,97	1064	445,06	419,98	82% / 100%